



HANGZHOU WEHERE MEDICAL TECHNOLOGY CO., LTD.

#701, Block C, Wanda Plaza, Gongshu, Hangzhou, 310015, China

Tel.: +86-571-86822701 Cell.: +8613905717534

E-mail: moneyxd@weheremed.com

www.weheremedical.com

WEHERE
MEDICAL
we here, we care



Anesthesia and Respiratory

+ Laryngeal Mask Airway	01
+ Endotracheal Tube	04
+ Breathing Filter & HME	09
+ Breathing Circuit Series	11
+ CPAP Tube	13
+ Masks	13
+ Incentive Spirometer	17
+ Aerosol Chamber	17
+ Airway Series	18
+ Manual Resuscitator	19
+ Infusion Pump	20

Surgical Drains

+ Drainage Tube	21
+ Suction Reservoir	24
+ Closed Wound Drainage System	26
+ Easy-Vac Suction Reservoir	29
+ Silicone Loop-Vessel Loops	29
+ Stainless Steel Trocar	30
+ Yankaure Handle	31

Urology & Gastroenterology

+ Stomach Tube	31
+ Urology Catheter	32
+ Suction Catheter	32
+ Silicone Gastrostomy Tube	33
+ Silicone Male External Catheter (Self Adhesive)	33

Other Products

+ Menstrual Cup	34
+ Menstrual Disc	35
+ Insulin Pen Needle	36
+ Blood Lancet (Round)	36
+ Safety Lancet	36

Disposable Silicone Laryngeal Mask Airway

For patients under anesthesia or pharmacological sedation and for patients requiring emergency artificial ventilation support during emergency and resuscitation to achieve upper airway patency

1. Made of medical grade silicone
2. Cuff is made of soft silicone and designed to conform to the contours of the throat, minimizing irritation and improving the seal with the throat.
3. Disposable & reusable both available, reusable type can be used up to 40 times by autoclave sterilization.



silicone check valve
for MRI compatible type



Colored pilot balloon
with size printed

Reusable Silicone Laryngeal Mask Airway



1. Bar/Aperture: With / Without
2. Cuff Color: Flesh Blue Transparent

Size	1.0	1.5	2.0	2.5	3.0	4.0	5.0	6.0
Single use Item No.	LM1004S1	LM1507S1	LM2010S1	LM2514S1	LM3020S1	LM4030S1	LM5040S1	LM6040S1
Single use Item No. (With Bar)	LM1004GS1	LM1507GS1	LM2010GS1	LM2514GS1	LM3020GS1	LM4030GS1	LM5040GS1	LM6040GS1
Reusable Item No.	LM1004S2	LM1507S2	LM2010S2	LM2514S2	LM3020S2	LM4030S2	LM5040S2	LM6040S2
Reusable Item No. (With Bar)	LM1004GS2	LM1507GS2	LM2010GS2	LM2514GS2	LM3020GS2	LM4030GS2	LM5040GS2	LM6040GS2
Patient Selection	0-5KG	5-10KG	10-20KG	20-30KG	30-50KG	50-70KG	70-100KG	70-100KG
Max. Inflation Volume	4ml	7ml	10ml	14ml	20ml	30ml	40ml	40ml

Disposable Reinforced Silicone Laryngeal Mask Airway

For patients under anesthesia or pharmacological sedation and for patients requiring emergency artificial ventilation support during emergency and resuscitation to achieve upper airway patency.

1. Made of medical grade silicone, non-toxic and non-irritating
2. High-strength reinforced wire inside, effectively prevent the tube body from flattening and kinking
3. Particularly suitable for upper face surgery, otolaryngology, ophthalmology and dental surgery. not interfere with the surgical field of view.
4. Disposable reinforced & reusable reinforced both available, reusable type can be used up to 40 times by autoclave sterilization.



Reusable Reinforced Silicone Laryngeal Mask Airway



Size	1.0	1.5	2.0	2.5	3.0	4.0	5.0	6.0
Single use Item No.	LM1004S1R	LM1507S1R	LM2010S1R	LM2514S1R	LM3020S1R	LM4030S1R	LM5040S1R	LM6040S1R
Single use Item No. (With Bar)	LM1004GS1R	LM1507GS1R	LM2010GS1R	LM2514GS1R	LM3020GS1R	LM4030GS1R	LM5040GS1R	LM6040GS1R
Reusable Item No.	LM1004S2R	LM1507S2R	LM2010S2R	LM2514S2R	LM3020S2R	LM4030S2R	LM5040S2R	LM6040S2R
Reusable Item No. (With Bar)	LM1004GS2R	LM1507GS2R	LM2010GS2R	LM2514GS2R	LM3020GS2R	LM4030GS2R	LM5040GS2R	LM6040GS2R
Patient Selection	0-5KG	5-10KG	10-20KG	20-30KG	30-50KG	50-70KG	70-100KG	70-100KG
Max. Inflation Volume	4ml	7ml	10ml	14ml	20ml	30ml	40ml	40ml
Remark	1. Bar/Aperture: With ✓ Without ✓ 2. Cuff Color: Flesh ✓ Blue ✓ Transparent ✓							

Disposable PVC Laryngeal Mask Airway

Disposable PVC Laryngeal Mask Airway is used for short-term artificial ventilation support for patients under general anesthesia or in need of emergency care, to ensure the patency of the upper airway.

1. Made of medical grade PVC, non-toxic and non-irritating.
2. Cuff is made of soft PVC and designed to conform to the contours of the throat.




Size	1.0	1.5	2.0	2.5	3.0	4.0	5.0
Item No.	LM1004P1	LM1507P1	LM2010P1	LM2514P1	LM3020P1	LM4030P1	LM5040P1
Reinforced Item No.	LM1004P1R	LM1507P1R	LM2010P1R	LM2514P1R	LM3020P1R	LM4030P1R	LM5040P1R
Patient Selection	0-5KG	5-10KG	10-20KG	20-30KG	30-50KG	50-70KG	70-100KG
Max. Inflation Volume	4ml	7ml	10ml	14ml	20ml	30ml	40ml
Remark	Only transparent color available, without bar.						

Disposable Preformed Silicone Laryngeal Mask Airway

1. 90 degrees flexible tubing is designed based on patient's laryngeal part of pharynx, which can be inserted comfortably, minimizing potential trauma and increasing sealing.

2. One piece design of cuff confirms the connection part softer.



Size	1.0	1.5	2.0	2.5	3.0	4.0	5.0
Item No.	LM1004SB2	LM1507SB2	LM2010SB2	LM2514SB2	LM3020SB2	LM4030SB2	LM5040SB2
Patient Selection	0-5KG	5-10KG	10-20KG	20-30KG	30-50KG	50-70KG	70-100KG
Max. Inflation Volume	4ml	7ml	10ml	14ml	20ml	30ml	40ml

Disposable/Reusable Double Lumen Reinforced Silicone Laryngeal Mask Airway

1. Upgrade double-lumen laryngeal mask airway.
2. Maintaining an effective airway.
3. Drainage lumen for suction the gastric juice.
4. Providing more airway security.



Size	1.0	1.5	2.0	2.5	3.0	4.0	5.0
Single use Item No.	LM1004S12	LM1507S12	LM2010S12	LM2514S12	LM3020S12	LM4030S12	LM5040S12
Reusable Item No.	LM1004S22	LM1507S22	LM2010S22	LM2514S22	LM3020S22	LM4030S22	LM5040S22
Patient Selection	<5KG	5-10KG	10-20KG	20-30KG	30-50KG	50-70KG	70-100KG
Max. Inflation Volume	4ml	7ml	10ml	14ml	20ml	30ml	40ml

Standard Endotracheal Tube (cuff/uncuffed)

1. Made of medical grade PVC, clear and smooth
2. High volum, low pressure cuff maintains good sealing
3. Murphy eye design, tube with X-ray opaque line


cuff

Item No.	ID(mm)	Item No.	ID(mm)
ET20SC	2.0	ET60SC	6.0
ET25SC	2.5	ET65SC	6.5
ET30SC	3.0	ET70SC	7.0
ET35SC	3.5	ET75SC	7.5
ET40SC	4.0	ET80SC	8.0
ET45SC	4.5	ET85SC	8.5
ET50SC	5.0	ET90SC	9.0
ET55SC	5.5	ET95SC	10.0



uncuffed

Item No.	ID(mm)	Item No.	ID(mm)
ET20S	2.0	ET65S	6.5
ET25S	2.5	ET70S	7.0
ET30S	3.0	ET75S	7.5
ET35S	3.5	ET80S	8.0
ET40S	4.0	ET85S	8.5
ET45S	4.5	ET90S	9.0
ET50S	5.0	ET95S	9.5
ET55S	5.5		
ET60S	6.0		




Reinforced Endotracheal Tube (cuff/uncuffed)

1. Made of medical grade PVC, clear and smooth.
2. High volume, low pressure cuff maintains good sealing.
3. Tube with spiral reinforced wire reducing risk of kinking.
4. Murphy eye design, tube with X-ray opaque line.


cuff

Item No.	ID(mm)	Item No.	ID(mm)
RET20PC	2.0	RET60PC	6.0
RET25PC	2.5	RET65PC	6.5
RET30PC	3.0	RET70PC	7.0
RET35PC	3.5	RET75PC	7.5
RET40PC	4.0	RET80PC	8.0
RET45PC	4.5	RET85PC	8.5
RET50PC	5.0	RET90PC	9.0
RET55PC	5.5	RET95PC	9.5



uncuffed

Item No.	ID(mm)	Item No.	ID(mm)
RET20P	2.0	RET60P	6.0
RET25P	2.5	RET65P	6.5
RET30P	3.0	RET70P	7.0
RET35P	3.5	RET75P	7.5
RET40P	4.0	RET80P	8.0
RET45P	4.5	RET85P	8.5
RET50P	5.0	RET90P	9.0
RET55P	5.5	RET95P	9.5




Endotracheal tube Oral Preformed (cuff/without cuff)

1. Oral preformed to reduce tube kinking during oral intubation.
2. Used for oral surgery, maxillofacial surgery.
3. Available with Magil Type.

cuff

Item No.	ID (mm)	Item No.	OD (mm)
OET20C	2.0	OET60C	6.0
OET25C	2.5	OET65C	6.5
OET30C	3.0	OET70C	7.0
OET35C	3.5	OET75C	7.5
OET40C	4.0	OET80C	8.0
OET45C	4.5	OET85C	8.5
OET50C	5.0	OET90C	9.0
OET55C	5.5	OET95C	9.5



uncuff

Item No.	ID (mm)	Item No.	OD (mm)
OET20	2.0	OET60	6.0
OET25	2.5	OET65	6.5
OET30	3.0	OET70	7.0
OET35	3.5	OET75	7.5
OET40	4.0	OET80	8.0
OET45	4.5	OET85	8.5
OET50	5.0	OET90	9.0
OET55	5.5	OET95	9.5




Endotracheal tube Nasal Preformed (cuff/without cuff)

1. Nasal preformed to reduce tube kinking during nasal intubation.
2. Used for nasal surgery, ophthalmic surgery, facial surgery.
3. Available with Magil Type.

cuff

Item No.	ID (mm)	Item No.	OD (mm)
NET20C	2.0	NET60C	6.0
NET25C	2.5	NET65C	6.5
NET30C	3.0	NET70C	7.0
NET35C	3.5	NET75C	7.5
NET40C	4.0	NET80C	8.0
NET45C	4.5	NET85C	8.5
NET50C	5.0	NET90C	9.0
NET55C	5.5	NET95C	9.5



uncuff

Item No.	ID (mm)	Item No.	OD (mm)
NET20	2.0	NET60	6.0
NET25	2.5	NET65	6.5
NET30	3.0	NET70	7.0
NET35	3.5	NET75	7.5
NET40	4.0	NET80	8.0
NET45	4.5	NET85	8.5
NET50	5.0	NET90	9.0
NET55	5.5	NET95	9.5




Endotracheal Tube with Suction Lumen

1. Made of medical grade PVC, clear and smooth.
2. High volum, low pressure cuff maintains good sealing.
3. Tube with suction lumen.
4. Murphy eye design, tube with X-ray opaque line.

Item No.	ID(mm)	Item No.	ID(mm)
SET060	6.0	SET080	8.0
SET065	6.5	SET085	8.5
SET070	7.0	SET090	9.0
SET075	7.5		

Item No.	ID(mm)	Item No.	ID(mm)
RSET060	6.0	RSET080	8.0
RSET065	6.5	RSET085	8.5
RSET070	7.0	RSET090	9.0
RSET075	7.5		



Double Lumen Endobronchial Tube

Used for thoracic surgery, bronchopulmonary manometry, administration of endobronchial anesthesia, and other common uses associated with sterile disposable bronchial intubation. Not for use in patients with primary bronchial stenosis or obstruction, where the primary bronchial anatomy would preclude proper tube placement.

1. Made of medical grade PVC.
2. Endobronchial Tube is included left-sided type and right-sided type.
3. It is composed of double-lumen bronchial intubation, multi-function connector, sputum suction tube and guide core.

Left-sided	Right-sided	Size(Fr)
EBT26L	EBT26R	26
EBT28L	EBT28R	28
EBT32L	EBT32R	32
EBT35L	EBT35R	35
EBT37L	EBT37R	37
EBT39L	EBT39R	39
EBT41L	EBT41R	41




Tracheostomy Tube (cuff/uncuffed)

1. Made of medical grade PVC, clear and smooth.
2. High volume, low pressure cuff maintains good sealing.
3. Tube with X-ray opaque line.
4. Round and smooth tip of obturator minimizes stimulation during intubation.

cuff

Item No.	ID(mm)	Item No.	ID(mm)
TT040PC	4.0	TT075PC	7.5
TT045PC	4.5	TT080PC	8.0
TT050PC	5.0	TT085PC	8.5
TT055PC	5.5	TT090PC	9.0
TT060PC	6.0	TT095PC	9.5
TT065PC	6.5	TT100PC	10.0
TT070PC	7.0		



uncuffed

Item No.	ID(mm)	Item No.	ID(mm)
TT025P	2.5	TT065P	6.5
TT030P	3.0	TT070P	7.0
TT035P	3.5	TT075P	7.5
TT040P	4.0	TT080P	8.0
TT045P	4.5	TT085P	8.5
TT050P	5.0	TT090P	9.0
TT055P	5.5	TT100P	10.0
TT060P	6.0		





Silicone Reinforced Endotracheal Tube


1. Both of the tube and cuff are made of medical grade silicone, better biocompatibility.
2. High volume, low pressure silicone cuff maintains good sealing.
3. Silicone tube with spiral reinforced wire reducing risk of kinking.
4. 100% medical grade silicone and soft tip, reduce the damage to the tracheal wall and airway mucosa.
5. Single-use and reusable 2 types available, both two types with stylet equipped.
6. Reusable type is autoclave sterilization only, can be reused up to 10 times.

Single-use type

Item No.	ID(mm)	OD(mm)
ET20PC	2.0	2.0
ET25PC	2.5	2.5
ET30PC	3.0	3.0
ET35PC	3.5	5.3
ET40PC	4.0	6.0
ET45PC	4.5	6.7
ET50PC	5.0	7.3
ET55PC	5.5	8.0
ET60PC	6.0	8.7
ET65PC	6.5	9.3
ET70PC	7.0	10.0
ET75PC	7.5	10.7
ET80PC	8.0	11.3
ET85PC	8.5	12.0
ET90PC	9.0	12.7
ET95PC	9.5	13.3

Reusable type


Item No.	ID(mm)	OD(mm)
ET20P	2.0	2.0
ET25P	2.5	2.5
ET30P	3.0	3.0
ET35P	3.5	5.3
ET40P	4.0	6.0
ET45P	4.5	6.7
ET50P	5.0	7.3
ET55P	5.5	8.0
ET60P	6.0	8.7
ET65P	6.5	9.3
ET70P	7.0	10.0
ET75P	7.5	10.7
ET80P	8.0	11.3
ET85P	8.5	12.0
ET90P	9.0	12.7
ET95P	9.5	13.3



Intubating Stylet

1. Made of malleable aluminum covered with plastic sheath.
2. Soft rounded tip for minimizing patient trauma during intubation
3. Flexible use for assistance of intubation


Item No.	Size(Fr)	Length	Applicable ETT
IS06	6	270mm	ID2.0-ID3.0
IS10	10	340mm	ID4.5-ID5.0
IS12	12	340mm	ID5.5-ID6.0
IS14	14	340mm	ID6.5-ID10.0



Endotracheal Tube Introducer

1. Made from semirigid non-toxic Polyethylene
2. Used for assistance intubation
3. Available with solid type and hollow type.
4. Hollow type can facilitate delivery of oxygen in emergency.

Size	Length (Solid type)	Length (Hollow type)
6ch	535mm	/
10ch	700mm	700mm
15ch	700mm	700mm



HMEF

HMEF is used for bacteria, particle filtration in breathing machine and anesthesia machine and to increase the gas moisture degree. The filter can effectively filter and stop the bacteria and virus between breathing system and breathing circuits, reduce the rate for lower respiratory tract infection, reduce the pain for patient and protect the equipment.



Item No.	Size	breathing volume(ml)	humidification efficiency	Resistance (pa)	Bacterial Filtration Efficiency(%)	Viral Filtration Efficiency(%)
CK-01	Adult	500	10mg/H2O@vt500ML	115@30L/min	99.99	99.99
CK-01-02	Adult	200-1500	/	100@30L/min	99.99	99.99
CK-02	Pediatric	500	10mg/H2O@vt502ML	215@30L/min	99.99	99.99
CK-02-01	Pediatric	150-300	24mg/H2O@Vt500ml	120@30L/min	99.99	99.99
CK-05	Neonate	70-150	22mg/H2O@Vt500ml	150@30L/min	99.99	99.99
CK-03-01	Adult	150-1500	32mg/H2O@Vt500ml	100@30L/min	99.99	99.99
CK-04-01	Adult	150-1500	19.5mg/500MLH2O	<0.3KPA(30L/MIN)	99.99	99.99
CK-41	Adult	150-1500	10mg/H2O@Vt500ml	13@30L/min	99.99	99.99
CK-001	Adult	100-1500	10mg/H2O@Vt501ml	100@30L/min	99.99	99.99

BVF & Spirometry Filter

This product can be equipped with pulmonary function machine to filter the bacterial from the patient. Available with mouthpiece and nose clip.

BVF

CKA-11-01



CKA-11



CKA-11-03



Spirometry Filter

CKC-34



CKC-39



Item No.	Size	ID (mm)	OD (mm)	Resistance (cm H ₂ O)	Bacterial Filtration Efficiency(%)	Viral Filtration Efficiency(%)
CKA-11-01	Adult	Standard	Standard	<1.5@60L/min	99.99	99.99
CKA-11	Adult	/	/	<1.3@60L/min	99.99	/
CKA-11-03	Adult	/	/	60pa	99.99	99.99
CKC-34	Adult	28.0/30.0/30.5	33.0	<1.3@60L/min	99.99	/
CKC-39	Adult	26.5	30.0	<1.3@60L/min	99.99	/
		29.5	33.5			

Tracheostomy HME

CK-24



CK-25



CK-26



CK-27



CK-28



Item No.	Size	Tidal volume	Humidity Output	flow resistance (pa)
CK-24	Adult	> 75ml	13.1mg/H ₂ O@vt500ML	15@30L/min
CK-25	Adult	> 75ml	8.9mg/H ₂ O@vt500ML	20@30L/min
CK-26	Adult	> 75ml	9.0mg/H ₂ O@vt500ML	50@30L/min
CK-27(Paper)	Adult	> 75ml	8.7mg/H ₂ O@vt500ML	10@30L/min
CK-27(Foam)	Adult	> 75ml	11.4mg/H ₂ O@vt500ML	15@30L/min
CK-28	Adult	> 75ml	7.8mg/H ₂ O@vt500ML	15@30L/min

Corrugated Circuit

Item No.	Size	Hose Length	Specification
BC1-L	Adult	1.2 / 1.6 / 1.8M	1. Standard connectors (15mm, 22mm). 2. EVA material, kinking resistant. 3. Normally 1.6m/1.8m, length can be customized
BC1-M	Pediatric	1.2 / 1.6 / 1.8M	
BC1-S	Infant	1.2 / 1.6 / 1.8M	



Extendable Circuit

Item No.	Size	Hose Length	Specification
BC2-L	Adult	1.2 / 1.6 / 1.8M	1. Standard connectors (15mm, 22mm). 2. EVA material, hoses are expandable, easy for transportation and use. 3. Normally 1.6m/1.8m, length can be customized.
BC2-M	Pediatric	1.2 / 1.6 / 1.8M	
BC2-S	Infant	1.2 / 1.6 / 1.8M	



Smoothbore Circuit

Item No.	Size	Hose Length	Specification
BC3-L	Adult	1.2 / 1.6 / 1.8M	1. Standard connectors (15mm, 22mm). 2. PVC material, smoothbore, kinking resistant, very high quality.
BC3-M	Pediatric	1.2 / 1.6 / 1.8M	
BC3-S	Infant	1.2 / 1.6 / 1.8M	



Coaxial Circuit

Item No.	Size	Hose Length	Specification
BC4-L	Adult	1.2 / 1.6 / 1.8M	1. Standard connectors (15mm, 22mm). 2. EVA material, Inside gas sampling line ends with swivel port. 3. Hose color: blue and purple for your choice
BC4-M	Pediatric	1.2 / 1.6 / 1.8M	
BC4-S	Infant	1.2 / 1.6 / 1.8M	



Midsplit Circuit

Item No.	Size	Hose Length	Specification
BC5-L	Adult	1.2 / 1.6 / 1.8M	1. Standard connectors (15mm, 22mm). 2. EVA material, single tube double chambers. 3. Normally 1.6m/1.8m, length can be customized.
BC5-M	Pediatric	1.2 / 1.6 / 1.8M	
BC5-S	Infant	1.2 / 1.6 / 1.8M	



Latex Breathing Bag



Latex-Free Breathing Bag




Silicone Breathing Bag




Breathing Circuit Set

Item No.	Type	Size
01	Breathing Tube	Adult - 1.5m,1.6m,1.8m, 2.0mAdult - 1.5m,1.6m,1.8m, 2.0m
02	Breathing Tube	Pediatric - 1.5m,1.6m,1.8m, 2.0m
03	Extra Lim Tube	1.0m,1.5m,1.8m, 2.0m
04	Breathing Bag	Adult - 3L,2L
05	Breathing Bag	Pediatric - 1L,0.5L
06	HMEF Filter	Adult
07	HMEF Filter	Pediatric
08	CO ₂ Monitoring Line	2m
09	Catheter Mount	
10	Catheter Mount Straight	
11	Y Connector with the CO ₂ Port	
12	Elbow Connector	
13	Luer Port & Cap 22/15mm	
14	Straight Connector 22mm	
15	Water Trap	



Anesthesia Circuit Set

Item No.	Type	Size
01	Breathing Tube	Adult - 1.5m,1.6m,1.8m, 2.0mAdult - 1.5m,1.6m,1.8m, 2.0m
02	Breathing Tube	Pediatric - 1.5m,1.6m,1.8m, 2.0m
03	Extra Lim Tube	1.0m,1.5m,1.8m, 2.0m
04	Breathing Bag	Adult - 3L,2L
05	Breathing Bag	Pediatric - 1L,0.5L
06	Anesthesia Mask	Adult - #4, #5, #6
07	Anesthesia Mask	Pediatric-#1, #2, #3
08	HMEF Filter	Adult
09	HMEF Filter	Pediatric
10	CO ₂ Monitoring Line	2m
11	Catheter Mount	
12	Catheter Mount Straight	
13	Y Connector with the CO ₂ Port	
14	Elbow Connector	
15	Luer Port & Cap 22/15mm	
16	Straight Connector 22mm	



Silicone Breathing Circuit - Reusable PC ⑪ PSF ⑬

1. Patent technology of integrated molding, high strength, not prone to falling off and being separated, good flexibility.
2. It can endure high temperature and high pressure disinfection, reusable.
3. Joints are made with injection molding with no gas leakage.
4. Mainly used for operation patients in anesthesia and oxygen patients after recovery.
5. Support and care for patient with severe postoperative respiratory.
6. Free options among water trap, Y type joint, L-shaped connector, masks, breathing bags, etc.




Item No.	Specification (ID)	Size
BCS-L	22mm	Adult
BCS-M	15mm	Pediatric
BCS-S	10mm	Neonate

Catheter Mount

1. Corrugated, smoothbore and extendable types are available for the hose.
2. Standard connectors (15F/22M, 22F), for adult and pediatric.

Item No.	Type	Size
CM-C01	Corrugated	Corrugated tube 150mm with Elbow 15F/22M and 22F Connector
CM-S01	Smoothbore	Smoothbore tube 150mm with Elbow 15F/22M and 22F Connector
CM-E01	Extendable	Extendable tube 150mm with Elbow 15F/22M and 22F Connector



CPAP Tube

1. Material: polypropylene (PP).
2. Connector ID: 22mm (Injection Molded).
3. Length: in general 1.8M, but can be customized (300mm- 3000mm).
4. Application: CPAP(Continuous Positive Airway Pressure), Sleep Apnea Tubing.



CPAP Mask

1. Made of medical grade PVC, soft and smooth.
2. For the treatment of obstructive sleep apnea by using continuous positive airway pressure.
3. Good Seal when fitted, providing comfort to patients.
4. PEEP Valve and Headgear can be optional equipped.
5. S, M, L three sizes available.

Item No.	Size (#)
CMP-S	S
CMP-M	M
CMP-L	L



CPR Mask

1. Made of medical grade PVC, soft and smooth.
2. Suitable for adult, children and infant.
3. Equip with a replaceable one-way valve with biological filter in one housing.
4. PVC gloves and alcohol pad can be optional equipped.
5. PP case and Nylon bag two packages available.



PVC Anesthesia Mask

1. The mask is made of medical grade PVC material.
2. The transparent mask shell is convenient for observing blood stain and natural respiration.
3. The contour of the air cushion of the mask fits the human face shape and provides good sealing effect.
4. The mask is disposable, reducing the risk of cross infection.

(Regular Type)



No Sterile Item No.	Sterile Item No.	Size (#)	Patient Size
AMPA0	AMPA0S	0	Neonate
AMPA1	AMPA1S	1	Infant
AMPA2	AMPA2S	2	Child Standard
AMPA3	AMPA3S	3	Small Adult
AMPA4	AMPA4S	4	Medium Adult
AMPA5	AMPA5S	5	Large Adult
AMPA6	AMPA6S	6	Extra Large Adult

PVC Anesthesia Mask

1. Made of medical grade PVC, clear and soft.
2. Seal tight against mouth.
3. Colored removable hook rings for easy size identification and fixation.

(Soft Cushion Type)



No Sterile Item No.	Sterile Item No.	Size (#)	Patient Size
AMPB0	AMPB0S	0	Neonate
AMPB1	AMPB1S	1	Infant
AMPB2	AMPB2S	2	Child Standard
AMPB3	AMPB3S	3	Small Adult
AMPB4	AMPB4S	4	Medium Adult
AMPB5	AMPB5S	5	Large Adult
AMPB6	AMPB6S	6	Extra Large Adult

PVC Anesthesia Mask (Easy Type)

1. Made of medical grade PP and TPE, clear and soft.
2. Soft anatomical cushion for good sealing and patient comfort.



No Sterile Item No.	Sterile Item No.	Size (#)	Patient Size
AMPC0	AMPC0S	0	Infant
AMPC1	AMPC1S	1	Infant
AMPC2	AMPC2S	2	Child Standard
AMPC3	AMPC3S	3	Small Adult
AMPC4	AMPC4S	4	Medium Adult
AMPC5	AMPC5S	5	Large Adult

Silicone Anesthesia Mask

Silicone material, good biocompatibility, silicone sealing edge makes the mask better fit the face.

1. The round design is suitable for infants, children's face shape.
2. High transparency, It is convenient to observe the patient.
3. Reusable.



One-piece Type



Split Type

Round Type




No Sterile Item No.	Sterile Item No.	Size (#)	OD
AMSC0	AMSC0S	0#	51.1mm
AMSC1	AMSC1S	1#	63.9mm
AMSC2	AMSC2S	2#	78.5mm

Oxygen Mask

1. Made of medical grade PVC.
2. Adjustable nose clip.
3. The length of tube could be customized.

No Sterile Item No.	Sterile Item No.	Size (#)
OMPXSS	OMPXS	XS
OMPSS	OMPS	S
OMPMS	OMPM	M
OMPLS	OMPL	L
OMPXLS	OMPXL	XL



Nebulizer Mask

1. To be used with oxygen tubing.
2. Nebulizer bottle volume: 6ml or 20ml.
3. With adjustable nose clip.
4. Standard connector.

Item No.	Size (#)
NMXSA	XS
NMSA	S
NMMA	M
NMLA	L
NMXLA	XL



Non-Rebreathing Mask

1. Made of clear, medical grade PVC.
2. With 1000/600/400/200ml reservoir bag.
3. Anti-crush tube, with elastic strap, with 2m oxygen tube.
4. Standard connector.
5. Latex free.

No Sterile Item No.	Sterile Item No.	Size (#)
OMNXSS	OMNXS	XS
OMNSS	OMNS	S
OMNMS	OMNM	M
OMNLS	OMNL	L
OMNXLS	OMNXL	XL



Nebulizer Kit With Mouthpiece

1. Made of clear, medical grade PVC and ABS.
2. Green/Clear to be selected.
3. Nebulizer bottle volume: 6ml or 20ml.
4. Smooth and feathered edge for patient comfort and reducing irritation points.

Item No.	Size (#)
NMXSB	XS
NMSB	S
NMMB	M
NMLB	L
NMXLB	XL



Tracheostomy Mask

1. Made of clear, medical grade PVC.
2. Swivel 360° adaptor with ID 18mm can connect with nebulizer.
3. Special snap design for easy removal.
4. Latex free.

No Sterile Item No.	Sterile Item No.	Size (#)
OMTXSS	OMTXS	XS
OMTSS	OMTS	S
OMTMS	OMTM	M
OMTLS	OMTL	L
OMTXLS	OMTXL	XL



Venturi Mask

1. Made of clear, medical grade PVC.
2. 7 density can be adjusted by two diluters.
3. Green: 24%,26%,28%,30%; White: 35%, 40%, 50%.
4. Adjustable nose clip, anti-crush tube.
5. With elastic strap.
6. With 2m oxygen tube.
7. Latex free.

No Sterile Item No.	Sterile Item No.	Size (#)
OMV2XSS	OMV2XS	XS
OMV2SS	OMV2S	S
OMV2MS	OMV2M	M
OMV2LS	OMV2L	L
OMV2XLS	OMV2XL	XL

With Two Diluters



Venturi Mask

1. Made of clear, medical grade PVC.
2. 7 diluters choice, each diluter for one oxygen density.
3. Adjustable nose clip.
4. Adaptor for high humidity.
5. Anti-crush tube.
6. Latex free.

Oxygen Density Setting	Color	Suggested Oxygen flow
24%	Blue	4LPM
28%	Yellow	4LPM
31%	White	6LPM
35%	Green	8LPM
40%	Pink	10LPM
50%	Orange	10LPM
60%	Red	12LPM

With Seven Diluters




No Sterile Item No.	Sterile Item No.	Size (#)
OMV6XSS	OMV6XS	XS
OMV6SS	OMV6S	S
OMV6MS	OMV6M	M
OMV6LS	OMV6L	L
OMV6XLS	OMV6XL	XL

Nasal Oxygen Cannula

1. Made of medical grade PVC, latex free.
2. 2m/2.1m oxygen tubing to be selected.
3. Soft prong, standard prong and flared prong three types available.

Item NO.	Size (#)	Patient Size
NOC01	S	Infant
NOC02	M	Pediatric
NOC03	L	Adult



ET/CO2 Sampling Cannula

1. Designed for oxygen supply and sampling of exhaled CO2 gas simultaneously.
2. Soft prongs minimize irritation.

Item NO.	Size (#)	Patient size
NOCL03	S	Infant
	M	Pediatric
	L	Adult



Incentive Spirometer

1. Deep Breathing Exerciser for comprehensive respiratory fitness - 3 chambers for differing inhalation rates from 600cc/sec to 1200 cc/ sec.
2. Keep your lungs and body healthy. The more breathing exercises performed properly, the more efficiently the lungs deliver oxygen to the entire body.
3. Lightweight to carry and clean, easy to put in a bag and carry it with you. You can use it at home, a simple exercise can keep your lungs healthy.

Item NO.	Size (#)
BSTH-600	600ml/sec
BSTH-900	900ml/sec
BSTH-1200	1200ml/sec
BSTH-700	700ml/sec
BSTH-1300	1300ml/sec
BSTH-1900	1900ml/sec



Aerosol Chamber

Aerosol Chamber is used to help deliver medication to patients who suffering from respiratory conditions such as chronic bronchitis. It is especially suitable for patients whose breathing and hand movements cannot coordinate, especially for young children and the elderly, and also for patients with pharyngitis or tonsillitis. It is small, easy to carry and has good effect, especially for children.

Item NO.	Patient Size
ICB-S	Infant
ICC-S	
ICB-M	Pediatric
ICC-M	
ICB-L	Adult
ICC-L	



Oropharyngeal Airway (Guedel Type)

The Guedel Airway is a non-tracheal catheter non-invasive ventilation pipe, which can prevent the tongue from falling back, quickly open the airway, and obtain effective ventilation through the middle passage.

1. The oropharyngeal airway is usually made of PE raw material.
2. It is easy to insert and can quickly establish airway.




Item No.	Size (#)	Length (mm)	Color Code
GA40S	000	40	Pink
GA50S	00	50	Blue
GA60S	0	60	Black
GA70S	1	70	White
GA80S	2	80	Green
GA90S	3	90	Yellow
GA100S	4	100	Red
GA110S	5	110	Light Blue
GA120S	6	120	Orange

Oropharyngeal Airway (Berman Type)

The Berman Airway is a non-tracheal catheter non-invasive ventilation pipe, which can prevent the tongue from falling back, quickly open the airway, and obtain effective ventilation from both sides of the passage.

1. The oropharyngeal airway is usually made of PE raw material.
2. It is easy to insert and can quickly establish airway.




Item No.	Size (#)	Length (mm)
BA40S	000	40
BA50S	00	50
BA60S	0	60
BA70S	1	70
BA80S	2	80
BA90S	3	90
BA100S	4	100
BA110S	5	110
BA120S	6	120

Nasopharyngeal Airway

Nasopharyngeal airway is suitable for patients with tight mandible, closed jaw and difficulty in inserting the oral airway.

1. Convenient and fast insertion.
2. Product specifications are complete. Single-use.



Item No. (green)	Item No. (blue)	Size(Fr)	ID(mm)
NAA25	NAB25	12	2.5
NAA30	NAB30	14	3.0
NAA35	NAB35	16	3.5
NAA40	NAB40	18	4.0
NAA45	NAB45	20	4.5
NAA50	NAB50	22	5.0
NAA55	NAB55	24	5.5
NAA60	NAB60	26	6.0
NAA65	NAB65	28	6.5
NAA70	NAB70	30	7.0
NAA75	NAB75	32	7.5
NAA80	NAB80	34	8.0
NAA85	NAB85	36	8.5
NAA90	NAB90	38	9.0
NAA95	NAB95	40	9.5

PVC Manual Resuscitator

It is applicable to the occasions of cardiopulmonary resuscitation and respiratory first aid at work. It is especially suitable for asphyxia, dyspnea or the need to increase the oxygen supply. It is easy to use, less pain, less complications, easy to carry, and can be immediately ventilated with or without oxygen source.

Conventional configuration: PVC resuscitation bag, air bag, oxygen tube, PVC anesthesia mask, instructions manual.

Optional accessories: oropharyngeal airway (edge type and central type), opener, PEEP valve.

Packing method: PP box, paper box, PE bag.

1. The product uses medical grade PVC. It is non-toxic and non-irritating.
2. The mask is made of PVC material, The air cushion contour of the mask fits the human face shape and provides good sealing effect.
3. The product specification is complete and is disposable



Item No.	Patient Size
PMRL	Adult
PMRM	Pediatric
PMRS	Infant

Silicone Manual Resuscitator



It is especially suitable for asphyxia, dyspnea or the need to increase the oxygen supply. It is easy to use, less complications, easy to carry, and can be immediately ventilated with or without oxygen source.

Conventional configuration: silicone resuscitation bag, air bag, oxygen tube, silicone anesthesia mask, instructions manual.

Optional 2 with 2 pieces: oropharyngeal airway (edge type and central type), opener, PEEP valve.

Packing method: PP box. Paper box, PE bag.

1. The product uses medical grade silicone, which is non-toxic, non-irritating and non-allergic.
2. The product specification is complete, and it is reused for 20 times.



Item No.	Patient Size
SMRL	Adult
SMRM	Pediatric
SMRS	Infant

Infusion Pump

1. For clinical pain-easing treatment, single use.
2. It uses the flexibility of the silicone reservoir to allow for a constant minimal drop control through the capillary catheter, guarantee a safe efficient continuous pain-easing infusion.
3. Used for continuous minimal infusion during clinical treatment.
4. Light-Shield infusion pump is available for photosensitive drugs.

CBI



Multirate, CBI



CBI+PCA



Multirate, CBI+PCA




Light-Shield infusion pump



Item No.	Type	Nominal Capacity (ml)	Basal Flow Rate (ml/hr)	PCA Bolus Volume	PCA Lockout Time
IP-C	CBI	60 / 100 / 150 / 200 / 275	2.0/3.0/4.0/5.0/10.0(normally2.0)	/	/
IP-MC	Multirate, CBI			0.5ml	15min
IP-P	CBI+PCA		2-4-6-8	/	/
IP-MP	Multirate, CBI+PCA		0-1-2-3-4-5-6-7 0-2-4-6-8-10-12-14	0.5ml	15min
Remark	Light-Shield infusion pump is available for photosensitive drugs.				

Silicone Round Fluted Drain (Hubless)



1. Made of imported medical grade silicone.
2. Utilize internal channels or fluted, to swiftly remove exudate.
3. Independent channels facilitate drainage and reduce risk of occlusion.
4. Made by one time, hubless with no glue.
5. Radio-opaque line through the length for X-ray visualization.
6. Available with or without "Three Face" stainless steel trocar.

Item No.	With Trocar	Size (Fr)	with connector
RFD10	RFD10T	10	XS
RFD12	RFD12T	12	S
RFD14	RFD14T	14	without
RFD15	RFD15T	15	without
RFD16	RFD16T	16	without
RFD18	RFD18T	18	without
RFD19	RFD19T	19	without
RFD22	RFD22T	22	without
RFD24	RFD24T	24	L
RFD28	/	28	L
RFD32	/	32	L

Silicone Round Fluted Drain



1. Made of medical grade silicone.
2. Round four channels shape allows for greater multi-channel drainage.
3. Channel drain offers important advantages over regular perforated drains by capillary effect.
4. Round Drains comes in a two-piece design.
5. Round Fluted part include silicone that is barium impregnated for X-ray detection.
6. Full fluted and 3/4 part fluted two types available.

Item No.	With Trocar	Fluted Drain (Fr)	Round Drain (Fr)	Type
RFD1001	RFD1001T	10	10	3/4 Part Fluted
RFD1501	RFD1501T	15	15	3/4 Part Fluted
RFD1901	RFD1901T	19	19	3/4 Part Fluted
RFD2401	RFD2401T	24	24	3/4 Part Fluted
FFD1002	FFD1002T	10	10	Full Fluted
RFD1502	RFD1502T	15	15	Full Fluted
RFD1902	RFD1902T	19	19	Full Fluted
RFD2402	RFD2402T	24	24	Full Fluted
Remark		Size and length can be customized.		


Silicone Flat Fluted Drain



1. It is made of medical grade silicone with good biocompatibility and is used for wound drainage.
2. The flat drainage tube consists of a silicone round tube and a white flat tube. The inner wall has a grooved structure to prevent the tube from blocking due to extrusion or twisting.
3. 3/4 fluted and full fluted two types available.
4. With radio-opaque function.
5. Stainless steel trocar can be used for surgical puncture, and suction reservoir can be used to collect drainage fluid for complete use.

Item No.	With Trocar	Flat Drain (mm)	Round Drain (Fr)	Type
FFD0701	FFD0701T	7	15	3/4 Part Fluted
FFD1001	FFD1001T	10	15	3/4 Part Fluted
FFD0702	FFD0702T	7	15	Full Fluted
FFD1002	FFD1002T	10	15	Full Fluted

Silicone Flat Perforated Drain



1. It is made of medical grade silicone with good biocompatibility.
2. The flat drainage tube consists of a transparent silicone round tube and a white flat tube. The inner wall has a toothed structure to prevent the tube from being blocked due to extrusion or twisting.
3. 3/4 perforated and full perforated two types available.
4. No connecting between flat tube and round tube, which reduces the patient's pain when removing the tube body.
5. With radio-opaque function.

Item No.	With Trocar	Flat Drain (mm)	Round Drain (Fr)	Type
FD0401	FD0401T	4	15	3/4 Part Perforated
FD0701	FD0701T	7	15	3/4 Part Perforated
FD1001	FD1001T	10	15	3/4 Part Perforated
FD1301	FD1301T	13	15	3/4 Part Perforated
FD0402	FD0402T	4	10	Full Perforated
FD0702	FD0702T	7	15	Full Perforated
FD1002	FD1002T	10	15	Full Perforated
FD1302	FD1302T	13	15	Full Perforated
Remark		Size and length can be customized.		



Silicone Round Perforated Drain



1. Made of imported medical grade silicone.
2. Radio-opaque line through the length for X-ray visualization.
3. Available with or without "Three Face" stainless steel trocar.

Item No.	With Trocar	Size (Fr)	with connector	Perforated segment
RPDA10	RPDA10T	10	XS	270±10 (mm)
RPDA12	RPDA12T	12	XS	
RPDA14	RPDA14T	14	S	
RPDA16	RPDA16T	16	/	
RPDA18	RPDA18T	18	/	
RPDA20	RPDA20T	20	/	
RPDA22	RPDA22T	22	L	
RPDA24	RPDA24T	24	L	
RPDA26	RPDA26T	26	/	
RPDA28	RPDA28T	28	/	
RPDA30	RPDA30T	30	L	200±10 (mm)
RPDA32	RPDA32T	32	L	
RPDB07	RPDB07T	7	2.4mm	
RPDB10	RPDB10T	10	3.5mm	
RPDB15	RPDB15T	15	/	
RPDB19	RPDB19T	19	/	
RPDB24	RPDB24T	24	/	




PVC Round Perforated Drain



1. Made of imported medical grade PVC.
2. Radio-opaque line through the length for X-ray visualization.
3. Available with or without "Three Face" stainless steel trocar.

Item No.	With Trocar	Size (Fr)
RPDP7	RPDP7T	7
RPDP8	RPDP8T	8
RPDP10	RPDP10T	10
RPDP12	RPDP12T	12
RPDP14	RPDP14T	14
RPDP15	RPDP15T	15
RPDP16	RPDP16T	16
RPDP18	RPDP18T	18
RPDP19	RPDP19T	19
RPDP20	RPDP20T	20
RPDP22	RPDP22T	22
RPDP24	RPDP24T	24



Silicone T Biliary Drain




1. Made of imported medical grade silicone.
2. T-shape silicone tubes for biliary drainage.
3. Total length is 35 cm and the horizontal part is 9.5 cm.



Item No.	Size (Fr)	T segment (mm)	Pipe section (mm)
TD08	8	95±5	325±20
TD09	9	97±5	
TD010	10	95±5	
TD012	12	97±5	
TD014	14	95±5	
TD015	15	97±5	
TD016	16	95±5	
TD018	18	97±5	
TD019	19	95±5	
TD020	20	97±5	
Remark	Size and length can be customized.		


Silicone Thoracic Drainage Tube



1. It is made of medical grade silicone .
2. X-ray line through the tube.
3. The front section of the catheter has smooth side holes to ensure smooth drainage.

Item No.	With Trocar	Size(Fr)	Item No.	With Trocar	Size(Fr)
TDT10	TDT10T	10	TDT26	TDT26T	26
TDT12	TDT12T	12	TDT28	TDT28T	28
TDT14	TDT14T	14	TDT30	TDT30T	30
TDT16	TDT16T	16	TDT32	TDT32T	32
TDT18	TDT18T	18	TDT34	TDT34T	34
TDT20	TDT20T	20	TDT36	TDT36T	36
TDT22	TDT22T	22	TDT38		38
TDT24	TDT24T	24	TDT40		40
Remark	Size and length can be customized.				

Silicone Penrose Drainage Tube



1. 100% medical grade silicone for superior biocompatibility.
2. Depend on gravitational forces, capillary action, and overflow for wounded fluid drainage.
3. Constant wall thickness over entire length.
4. X-ray opaque line throughout the length.
5. Clear connecting tube provides inspection and drainage visibility.

Item No.	Size (mm)	Drain type	Length
RTF0630	6	Flat	30 cm
RTF0830	8		
RTF1030	10		
RTF1230	12		
RTF1430	14		
RTF1630	16		
RTF2530	25	Flat	50 cm
RTF4030	40		
RTF0650	6		
RTF0850	8		
RTF1050	10		
RTF1250	12		
RTF1450	14		
RTF1650	16		
RTF2550	25		
RTF4050	40		

Round Drainage Tube



The material is translucent silicone material for medical tubing and profiles. Our tubing material offers very good mechanical properties, moderate solvent resistance and excellent temperature tolerance.

1. 100% medical grade silicone.
2. Clear tube provides inspection and drainage visibility.
3. Non-sterile silicone tubings.

Item No.	Size (mm)	CH/FR	Item No.	Size (mm)	CH/FR	Item No.	Size (mm)	CH/FR
SRT0701	1.6X2.4	CH7.2	SRT2103	5X7	CH21	SRT3601	7X12	CH36
SRT0901	1 X3	CH9	SRT2401	5X8	CH24	SRT3602	8X12	CH36
SRT0902	1.5X3	CH9	SRT2402	6X8	CH24	SRT3901	7X13	CH39
SRT1201	2X4	CH12	SRT2701	5X9	CH27	SRT3902	9X13	CH39
SRT1501	3X5	CH15	SRT2702	6X9	CH27	SRT4201	8X14	CH42
SRT1801	3X6	CH18	SRT3001	6X10	CH30	SRT4202	10X14	CH42
SRT1802	4X6	CH18	SRT3002	7X10	CH30	SRT4501	10X15	CH45
SRT2101	3X7	CH21	SRT3301	7X11	CH33	SRT4801	10X16	CH48
SRT2102	4X7	CH21	SRT3302	8X11	CH33	SRT5401	10X18	CH54
						SRT5701	10X19	CH55

Suction Reservoir



1. Made of imported medical grade silicone.
2. Provided with internal anti-reflux valve.
3. Bulb silicone walls provide easy activation of suction.
4. Luer cap design, more safe and convenient to process when using.
5. Regular system is to be used together with silicone flat drains, but we also recommend innovative solutions which are combined with our other silicone drains.

Item No.	Size(ml)	Graduation(ml)
SR1100	100	25/50/75/100
SR1150	150	25/50/75/100/125/150
SR1200	200	50/100/150/200
SR1400	400	100/200/300/400

Suction Reservoir Bottom Outlet

1. Silicone bulbs are suitable for plastic surgery, neurosurgery, orthopaedics, OB/GYN surgery and laparoscopic surgery.
2. Bottom outlet design are available for size 100ml and 200ml silicone reservoir.

A Ruhr hat	Cap stopper	Size(ml)	Graduation(ml)
SR4100A	SR5100B	100	25/50/75/100
SR4150A	SR5150B	150	50/100/150
SR4200A	SR5200B	200	50/100/150/200
SR4400A	SR5400B	400	50/100/150/200



Drainage Bag

1. The drainage bag is designed to be connected with our silicone reservoir to collect, store and record the fluid after the surgeries.
2. Be equipped with internal one-way valve to prevent reflux.
3. Be equipped with an extension PVC tube and a luer lock to be connected with our silicone reservoir.
4. Can be replaced individually.

Item No.	Size(ml)
DB11000A	1000
DB11000B	1000



Closed Wound Drainage System (With Trocar)



1. This product includes silicone reservoir, silicone flat drains or other silicone drains and stainless steel trocar.
2. Can use separately or using trocar for insertion and combining with silicone reservoir for fluids suction and collection.
3. Provided with internal anti-reflux valve to prevent backflow of fluid to the patient.
4. Silicone bulbs are suitable for plastic surgery, neurosurgery, orthopaedics, OB/GYN surgery and laparoscopic surgery.

Item No.	Size(ml)	Type
SDS100FD1 T	100	3/4 Part Fluted
SDS150FD1 T	150	3/4 Part Fluted
SDS200FD1 T	200	3/4 Part Fluted
SDS400FD1 T	400	3/4 Part Fluted
SDS100FD2T	100	Full Fluted
SDS200FD2T	200	Full Fluted

Wound Drainage system (with Drainage Bag)



1. This product includes silicone reservoir, silicone round drains or other silicone drains, stainless steel trocar and drainage bag.
2. Silicone bulbs are suitable for plastic surgery, neurosurgery, orthopaedics, OB/GYN surgery and laparoscopic surgery.
3. The drainage bag is designed to be connected with our silicone bulb reservoir to collect, store and record the fluid after the surgeries.
4. The drainage bag with the capacity of 1000ml, is equipped with internal one-way valve to prevent reflux.

Item No.	Size(ml)	Type
SDS100FD1 T	100	3/4 Part Fluted
SDS150FD1 T	150	3/4 Part Fluted
SDS200FD1 T	200	3/4 Part Fluted
SDS400FD1 T	400	3/4 Part Fluted
SDS100FD2T	100	Full Fluted
SDS200FD2T	200	Full Fluted

Suction Reservoir (Bottom outlet with drainage bag)



1. This product includes silicone reservoir (bottom outlet), silicone round drains or other silicone drains, stainless steel trocar and drainage bag.
2. Silicone bulbs are suitable for plastic surgery, neurosurgery, orthopaedics, OB/GYN surgery and laparoscopic surgery.
3. The drainage bag is designed to be connected with our silicone bulb reservoir to collect, store and record the fluid after the surgeries.
4. The drainage bag with the capacity of 1000ml, is equipped with internal one-way valve to prevent reflux.

Item No.	Size(ml)
SR4100A	100
SR4150A	150
SR4200A	200
SR4400A	400

Closed Wound Drainage System

1. This product includes silicone reservoir (bottom outlet), silicone round drains or other silicone drains, stainless steel trocar and drainage bag.
2. Silicone bulbs are suitable for plastic surgery, neurosurgery, orthopaedics, OB/GYN surgery and laparoscopic surgery.
3. The drainage bag is designed to be connected with our silicone bulb reservoir to collect, store and record the fluid after the surgeries.
4. The drainage bag with the capacity of 1000ml, is equipped with internal one-way valve to prevent reflux.

Item No.	Size(ml)	Type
SDS100RFDT	100	With Trocar
SDS150RFDT	150	With Trocar
SDS200RFDT	200	With Trocar
SDS400RFDT	400	With Trocar



Closed Wound Drainage System (Spring with silicone drain)

1. This product including 3-spring evacuator, silicone connecting tubing with
2. Y connector, silicone round perforated drains with trocar.
3. The design of 3-spring evacuator minimizes the Collapse-Resistance.

Item No.	Size
PDS200RFDT	200ml+15Fr
	200ml+19Fr
PDS400RFDT	400ml+15Fr
	400ml+19Fr



Closed Wound Drainage System (Spring with PVC drain)

1. This product including 3-spring evacuator, PVC connecting tubing with
2. Y connector, PVC round perforated drains with trocar.
3. The design of 3-spring evacuator minimizes the Collapse-Resistance.

Item No.	Size(ml)	Item No.	Size(ml)
PDS200RPDPT	200ml+7Fr	PDS400RPDPT	400ml+7Fr
	200ml+10Fr		400ml+10Fr
	200ml+12Fr		400ml+12Fr
	200ml+14Fr		400ml+14Fr
	200ml+16Fr		400ml+16Fr
	200ml+18Fr		400ml+18Fr



Closed Wound Drainage System (Hollow)

1. Aseptic PE material.
2. Hollow non-spring wound drainage reservoir with PVC drainage tube and trocar, bed sheet clamp and robot clamp.

Item No.	Size(ml)	Item No.	Size(ml)
PDH 25ML	25ml+10Fr	PDH 100ML	100ml+10Fr
	25ml+12Fr		100ml+12Fr
	25ml+14Fr		100ml+14Fr
	25ml+16Fr		100ml+16Fr
	25ml+18Fr		100ml+18Fr
Item No.	Size(ml)	Item No.	Size(ml)
PDH 400ML	400ml+10Fr	PDH 800ML	800ml+10Fr
	400ml+12Fr		800ml+12Fr
	400ml+14Fr		800ml+14Fr
	400ml+16Fr		800ml+16Fr
	400ml+18Fr		800ml+18Fr



Straight Connector

1. Fits 3/8 inch - 1/2 inch x 5/32 inch - 5/16 inch ID Tubing (0.375 inch - 0.5 inch x 0.156 inch - 0.312 inch; 9.5mm - 12.7mm x 4mm - 7.9mm)
2. Made of PA material.
3. Connector can be customized upon request.



Item No.

WDA01

Adaptor for Drains

1. Made of PP material.

Item No.	Size	Application	Item No.	Size	Application
WDA02	XS	10/12Fr	WDA05	1#	7-8 Fr
WDA03	S	14Fr	WDA06	2#	10-12 Fr
WDA04	L	22/24/28Fr	WDA07	3#	15-19 Fr
			WDA08	4#	24-28Fr




Easy-Vac Suction Reservoir



Item No.	Size(ml)
EV200	200
EV400	400

1. Closed System Design reduces the chance of cross contamination.
2. Compact and Self Contained System allows complete patient ambulation.
3. Pre-compressed State of Package saves the shipping cost and storage space.
4. Y-connector for Bilateral Drainage.
5. Attached Stepped Adapters allows connecting to other sizes on demand which facilitates easy connection to drainage system.
6. Internal Anti-reflux Valve is incorporated to prevent reverse flow during emptying and reactivation.
7. Reclosable Pouring Spout permits convenient emptying of the reservoir unit and system reactivated when required.
8. Latex-free Construction eliminates latex-induced allergic reactions. Reservoir is available in 2 sizes (200 & 400 ml).

Silicone Loops



1. Four bright colors for easy identification of vessels, veins, nerves and tendons throughout the procedur.
2. Made from durable, fluid-repellent, radiopaque, medical-grade silicone for supreme durability.
3. Neatly wrapped to prevent tangling.

Item No.	Description	Item No.	Description
VPA01	Blue Micro	VPA03	Blue Maxi
	Red Micro		Red Maxi
	White Micro		White Maxi
	Yellow Micro		Yellow Maxi
VPA02	Blue Mini	VPA04	Blue SuperMaxi
	Red Mini		Red SuperMaxi
	White Mini		White SuperMaxi
	Yellow Mini		Yellow SuperMaxi

Stainless Steel Trocar

1. Made of stainless steel 304 material.

6 inch Trocar

Item No.	Size	Item No.	Size
T6-10	10Fr	T6-19	19Fr
T6-12	12Fr	T6-20	20Fr
T6-14	14Fr	T6-22	22Fr
T6-15	15Fr	T6-24	24Fr
T6-16	16Fr		


8 inch Trocar

Item No.	Size	Item No.	Size
T8-10	10Fr	T8-18	18Fr
T8-12	12Fr	T8-19	19Fr
T8-14	14Fr	T8-20	20Fr
T8-15	15Fr	T8-22	22Fr
T8-16	16Fr	T8-24	24Fr



6/8 inch Trocar

Item No.	Size	Item No.	Size
TH8-07	7Fr	TH8-15	15Fr
TH8-10	10Fr	TH8-16	16Fr
TH8-12	12Fr	TH8-18	18Fr
TH8-14	14Fr	TH8-20	20Fr



Hook Type Trocar


Bendable Trocar

Item No.	Size	Item No.	Size
TB6-10	10Fr	TB6-19	19Fr
TB6-15	15Fr	TB6-24	24Fr



Double Bendable Trocar

Item No.	Size	Item No.	Size
TDB8-10	10Fr	TDB8-15	15Fr
TDB8-12	12Fr	TDB8-19	19Fr
TDB8-14	14Fr	TDB8-24	24Fr



Yankauer Handle

1. Made of transparent material for good visualization during suction.
2. The forming angle of handle is designed to make it more convenient and humanized in clinic use.

Item No.	Size
YH01	StandardTip, with control Vent
YH02	StandardTip, without control Vent
YH03	CrownTip, with control Vent
YH04	CrownTip, without control Vent
YH05	RoundTip, with control Vent
YH06	RoundTip, without control Vent



Suction Tubing With Yankauer Handle

1. Made of transparent material.
2. Can be combined with different types of yankauer handle and connecting tube to suit various requirements.


Item No.	Size
YH01SV18	StandardTip, with control Vent + 3/16", 1.8M
YH01SW36	StandardTip, without control Vent + 3/16", 3.6M
YH02C18	CrownTip, with control Vent + 1/4", 1.8M
YH02C36	CrownTip, without control Vent + 1/4", 3.6M
YH03R18	RoundTip, with control Vent + 9/32", 1.8M
YH03R36	RoundTip, without control Vent + 9/32", 3.6M



Silicone Stomach Tube

1. Made of 100 % medical grade silicone.
2. Provision of four eye-holes for efficient aspiration
3. With X-ray line, 120CM.


Normal Tip		Ryle Tip	
Item No.	Size (Fr)	Item No.	Size (Fr)
ST06S	6	ST08SR	8
ST08S	8	ST10SR	10
ST10S	10	ST12SR	12
ST12S	12	ST14SR	14
ST14S	14	ST16SR	16
ST16S	16	ST18SR	18
ST18S	18	ST20SR	20
ST20S	20	ST22SR	22
ST22S	22	ST24SR	24
ST24S	24		
ST36S	36		
ST39S	39		



PVC Stomach Tube

1. Made of 100 % medical grade PVC.
2. Provision of four eye-holes for efficient aspiration
3. With X-ray line, 120CM.

Normal Tip			
Item No.	Size (Fr)	Item No.	Size (Fr)
ST06P	6	ST16P	16
ST08P	8	ST18P	18
ST10P	10	ST20P	20
ST12P	12	ST22P	22
ST14P	14	ST24P	24



Silicone Foley Catheter

1. Made of 100% medical grade silicone.
2. Symmetrical balloon expands equally so as to work safely and efficiently.
3. Radio-opaque line through the length for X-ray visualization.

2-Way

Item No.	Size (Fr)	Balloon Volume (ml)
FC06S2	6	3
FC08S2	8	3/5
FC10S2	10	3/5
FC12S2	12	5/10/15/30
FC14S2	14	5/10/15/30
FC16S2	16	5/10/15/30
FC18S2	18	5/10/15/30
FC20S2	20	5/10/15/30
FC22S2	22	5/10/15/30
FC24S2	24	5/10/15/30

3-Way

Item No.	Size (Fr)	Balloon Volume (ml)
FC16S3	16	30
FC18S3	18	30
FC20S3	20	30
FC22S3	22	30
FC24S3	24	30
FC26S3	26	30



PVC Nelaton Catheter

1. Tiemann tip available.
2. Hydrophilic coated type available.
3. Both female and male types available.

Item No.	Size (Fr)	Item No.	Size (Fr)
NC06	6	NC16	16
NC08	8	NC18	18
NC10	10	NC20	20
NC12	12	NC22	22
NC14	14	NC24	24



Suction Catheter

Item No.	Size	Item No.	Size
SC06	6Fr	SC14	14Fr
SC08	8Fr	SC16	16Fr
SC10	10Fr	SC18	18Fr
SC12	12Fr	SC20	20Fr




Closed Suction Catheter

1. Color-coded control valve , easy to identify sizes.
2. Designed to clear bacterial build up and reduce VAP pathogenesis.
3. Used with Endotracheal Tube or Tracheotomy Tube.
4. 24 hours and 72 hours two types are available.

Normal Tip Type

Item No.	Size (Fr)	Size of Endotracheal Tube/ Tracheotomy Tube: I.D.(mm)
CSC06	6	4.0-4.5
CSC08	8	5.0-6.0
CSC10	10	6.5-7.0
CSC12	12	7.5-8.0
CSC14	14	8.5-9.0
CSC16	16	9.5-10.0



Silicone Gastrostomy Tube

1. Suitable for gastrostomy.
2. Made of medical silicone, have good biocompatibility.
3. The tube have large lumen can effectively reduce the tube occlusion.
4. Have Radio-opaque line for detection of correct placement.
5. Short catheter design help the balloon close to the stomach wall, have good elasticity and flexibility, can reduce stomach trauma.
6. The multi-function connector have Feeding Port and Medication Port provide a variety of connecting use more easily and quickly.
7. Color coding for size identification


REF	Size	Balloon Volume	Color Code	OD	Length
SG0112	12Fr	3-5ml	White	4.0mm	235mm
SG0114	14Fr	3-5ml	Green	4.7mm	235mm
SG0116	16Fr	5-20ml	Orange	5.3mm	235mm
SG0118	18Fr	5-20ml	Red	6.0mm	235mm
SG0120	20Fr	5-20ml	Yellow	6.7mm	235mm
SG0122	22Fr	10-20ml	Purple	7.3mm	235mm
SG0124	24Fr	10-20ml	Blue	8.0mm	235mm



Silicone Male External Catheter (Self Adhesive)

1. All silicone construction provides excellent breath-ability. Available size selection helps ensure a secure fit. LATEX FREE Safe, Strongly Skin Friendly Adhesive Bonding (Good Initial cohesiveness, and good holding power time)
2. Medical Grade Solution Acrylic Pressure Sensitive Adhesive was applied to self adhesive area which make patient more confidence when they using the products with Urine bag. It could be kept 24 hours without any dropping problems. The lab test show the peeling strength could be above 10N/INCH

Styles	Sizes				
Pop On	24mm	28mm	31mm	35mm	40mm
Standard	24mm	28mm	31mm	35mm	40mm
Wide Band	24mm	28mm	31mm	35mm	40mm




Menstrual Cup

The menstrual cup is reusable and made of medical grade liquid silicone. It is designed to be folded and inserted into the vagina to collect menstrual blood. Once in, it will open, lower than the tampon.

1. Made of medical grade liquid silicone.
2. Size: T, S, L three sizes available.
3. Color: Transparent, Pink, Purple, Green and Rose Red.
4. Luke Cup and Vivi Cup two types available.
5. Customizable packaging according to customer requirements.


Luke Cup

Item No.	size	Total Length (B)	Total Length (D)	Capacity to the holes
MC1601	T	43.1±0.5 mm	38.3±0.5 mm	17ml
MC1602	S	48.2±0.5 mm	42.5±0.5 mm	24ml
MC1603	L	53.8±0.5 mm	45.2±0.5 mm	30ml



Vivi Cup

Item No.	Size	Total Length (B)	Total Length (D)	Capacity to the holes
MC1901	T	39.7±1.0 mm	38.3±0.5 mm	16ml
MC1902	S	42.7±1.0 mm	42.0±0.5 mm	21ml
MC1903	L	47.2±1.0 mm	46.0±1.0 mm	26ml



Optional Accessories

Foldable Container

1. Made of 100% food grade silicone, microwave compatible.
2. 8.1× 8.7 cm (unfolded), 8.1x 2.4 cm (olded).
3. Green and orange two colors available.

Manual

1. 100% recycled paper, A4 size.
2. Folded manual size is 5.5 x 7.3 cm

Bamboo Brush

1. Bamboo handle with nylon brush.
2. 8x 1.1x 1.7 cm.

Carry Bag

1. Microfiber and 100% cotton two types available.
2. 10 x 13 cm.

Manual



Bamboo Brush



Foldable Container



Carry Bag



Color Box -A



Color Box -B



Color Box -C



Color Box -D




Color Box -E

Menstrual Disc

1. Each disc is made of medical grade silicone.
2. Greater comfort, greater capacity, no additives, no smell.
3. Can be used for up to 12 consecutive hours.


Menstrual Disc 001




Menstrual Disc 002



Menstrual Disc 003



Menstrual Disc 004



Menstrual Disc 005



Menstrual Disc 006



Insulin Pen Needle

1. Trip-bevel needle tip and thin wall technology reduced the injection resistance to a large extent, minimize the injury to the skin. Pain-free experience during injection, special needle tip design can effectively avoid the needle bending or breaking during usage.

2. Ultra fine needle (33G) with shortest length(4mm) reduce the muscle layer risk, preforming a more safeand eaiser injection.

3. The ultra-thin wall needle tube design improves the flow, ensures smooth injection, and reduces injection time and force.

4. The double layer silicone coating technology of the needle tube makes the injection more smooth and less painful.

Item No.	Metric System	Metric System	Color
IPNO13304	0.20x4mm	33Gx4mm	Green
IPNO13204	0.23x4mm	32Gx4mm	Green
IPNO13205	0.23x5mm	32Gx5mm	Purple
IPNO13206	0.23x6mm	32Gx6mm	Blue
IPNO13208	0.23x8mm	32Gx8mm	Red
IPNO13104	0.25x4mm	31Gx4mm	Green
IPNO13105	0.25x5mm	31Gx5mm	Purple
IPNO13106	0.25x6mm	31Gx6mm	Blue
IPNO13108	0.25x8mm	31Gx8mm	Red
IPNO13006	0.30x8mm	30Gx8mm	Yellow
IPNO12912	0.33x12mm	29Gx12mm	Purple



Blood Lancet (Round)

Twist off lancets are manufactured with the goal of providing a simple, safe and easy method to achieve capillary blood sampling. The lancets can be used with most currently available lancing devices. Twist lancets are available in gauge sizes from 21G to 33G

Item No.	Size	Color	Package Type
BL0221	21G	Yellow	25pcs/box, 5000pcs/carton 50pcs/box 10000pcs/carton 100pcs/box 20000pcs/carton
BL0223	23G	Blue	
BL0226	26G	Purple	
BL0228	28G	Blue	
BL0230	30G	Purple	
BL0232	32G	Yellow	
BL0233	33G	Green	



Safety Lancet

1. Directly touch the sampling site and press, no need to use lancing device.

2. Two springs design and tri-beveled needle tip ensure high speed penetration and minimize the pain, which makes blood sampling easy and convenient.

3. Independently development, patented technology. Self-destruction structure design let the medical staff and patients feel more safe and reliable.

Item No.	Size	Color	Packing
SLO0321	21G/1.8mm	Orange	50boxes/ 100boxes/ 1carton
SLO0323	23G/1.8mm	Grey	
SLO0326	26G/1.8mm	Yellow	
SLO0328	28G/1.8mm	Purple	
SLO0330	30G/1.8mm	Teal	

