



HANGZHOU WEHERE MEDICAL TECHNOLOGY CO., LTD.

#701, Block C, Wanda Plaza, Gongshu, Hangzhou, 310015, China

Tel.: +86-571-86822701 Cell.: +8613905717534

E-mail: moneyxd@weheremed.com

www.weheremedical.com

WEHERE
MEDICAL
we here, we care



CHAPTER 1

RESPIRATORY

EXCELLENCE QUALITY INTEGRITY PROFESSIONAL

⌘ PRODUCTS

RESPIRATORY PRODUCTS

Non-Rebreating Oxygen Mask

Features

- 1.Mask is made of medical grade PVC soft & comfortable
- 2.Adjustable nose clip
- 3.Elastic Strap,latex-free is optional
- 4.Reservoir bag of 600ml / 800ml /1000ml
- 5.7 fts Oxygen tube
- 6.Available with transparent / green color
- 7.Customization is acceptable



Nebulizer Mask

Features

- 1.Mask is made of medical grade PVC soft & comfortable
- 2.Adjustable nose clip
- 3.Elastic Strap, latex-free is optional
- 4.Nebulizer bottle volumes of 6ml, 8ml, 20ml
- 5.7 fts oxygen tube
- 6.Available with transparent / green color
- 7.Customization is acceptable



Adjustable Venturi Mask 1

Features

- 1.Mask is made of medical grade PVC soft & comfortable
- 2.Adjustable nose clip
- 3.Elastic Strap, latex-free is optional
- 4.Oxygen concentration range: 24-50%
- 5.7 fts Oxygen tube
- 6.Available with transparent / green color
- 7.Customization is acceptable



Adjustable Venturi Mask 2

Features

- 1.Mask is made of medical grade PVC soft & comfortable
- 2.Adjustable nose clip
- 3.Elastic Strap,latex-free is optional!
- 4.Color-coded diluters for easy identification of oxygen concentration
- 5.Oxygen concentration range: 24-60%
- 6.7 fts Oxygen tube
- 7.Customization is acceptable



Nasal Oxygen Cannula

Features

- 1.Made of medical grade PVC soft & comfortable
- 2.Available with injection prongs and soft prongs
- 3.Anti-crush oxygen delivery tube
- 4.Available with transparent / green color
- 5.Tube length: 7 fts
- 6.Customization is acceptable



CO2 Monitoring Cannula

Features

- 1.Simultaneous CO2 monitoring and oxygendelivery
- 2.CO2 sampling through nasal prongs foraccurate monitoring
- 3.Tubing length: 7 fts
- 4.Available with transparent / green color
- 5.Customization is acceptable



Oxygen Mask

Features

- Material: Medical grade PVC
Components: With adjustable nose clip, elastic strap, 7 fts oxygen tube
Size: S, M, L, XL
Color: Green, white transparent



Breathing Exerciser / Spirometer

Features

- Size:500-5000CC
Type:Single,Two,Three Ball
Packing:1pc/box
Material:ABS
Color:Custom



Heat And Moisture Exchange Filters (HMEF)

Features

- 1. Effectively filter and stop 99.99% of bacteria and virus between breathing system and breathing circuits.
 - 2. Adjust the humidity of the air patients take in.
 - 3. Exchange the thermal energy of air breath in and out.
- Material: ABS, filter paper
Size: Adult / Pediatric



CO2 Sampling Line With Filter

Features


- 1.Made of medical grade PVC material;
- 2. Tube length 3 meters, customization is acceptable;
- 3. Sampling lines for measurement of end-tidal CO2;
- 4. The filter protects monitor from moisture and contamination;
- 5. Male & Female Luer locks are compatible with most of the monitoring devices on the market.



Breathing Bag

Features

- The reservoir bag is used for anesthesia and breathing circuits. If latex material is allergic to patient, then latex-free artificial breathing apparatus should be used.
- 1. Material: Latex/ Latex-free
 - 2. A single patient use, can reduce cross contamination
 - 3. Bag volume: 3L/2L/1L/0.5L



Central Venous Catheter Kit

Central Venous Catheter

Kits Include:

- | | |
|-----------------------------|-----------------------|
| 1. Central Venous Catheter | 6. Introducer Needle |
| 2. Guidewire | 7. Introducer Syringe |
| 3. Vessel Dilator | 8. Injection Needle |
| 4.Clamp | 9. Injection Cap |
| 5. Fastener: Catheter Clamp | |



CHAPTER 2

ANESTHESIA

EXCELLENCE QUALITY INTEGRITY PROFESSIONAL

⌘ PRODUCTS

ANESTHESIA PRODUCTS

Anesthesia Breathing System

Features

Application:

Connect the equipment, transfer the air to body

Size: Adult, pediatrics

Type: common, extension tube

Construction with L connector, three-way connector

Attention: Do not give heavy press on it when store



Anesthesia Mask

Features

Material: TPE & PP

Advantages: No smell, no need assemble
& can be used only by one hand.

Size: 0/1/2/3/4/5/6



Closed Suction System

Features

1. Made of medical grade material, easy operation
2. Roundness & Smooth of the patients' end of catheter, minimize the injury of respiratory mucosa
3. Transparent 5-way design, oxygen delivery & suction performed simultaneously
4. The cover highly sealed, effectively prevent the cross-infection
5. Built-in spring of vacuum control, finger-on for suction, finger-off for closing
6. 7-tone date indication to remind the replacement
7. Size available from Fr 6 - Fr 16
8. Type: 24hrs, 72hrs



Tracheostomy Tube

Features

1. Made from medical grade PVC.
2. Supplied sterile in individual peelable blister pack.
3. For single use only.
4. Size: 3.0 - 9.5 mm
5. Type: cuff / plain



Endotracheal Tube

Features

- 1.Made of non-toxic medical-grade PVC, transparent,soft and smooth.
- 2.With high volume low press cuff.
- 3.Radio opaque line through the length for x-ray visualization.
- 4.Availabel with Murphy eye only.
- 5.Size:2.0mm-10.0mm.
- 6.Type: cuff/plain.



Reinforced Endotracheal Tube

Features


- 1.Made of non-toxic medical-grade PVC, transparent, soft and smooth
- 2.With high volume low press cuff.
- 3.Spiral reinforcement minimizes crushing or kinking.
- 4.Conform to any patient posture,especially to operation of decubitus.
- 5.Size:2.5m-10.0mm.
- 6.Type:cuff/plain



Endotracheal Tube With Aspiration Lumen (standard and reinforced type)

Features


- 1.Effective suctioning of accumulated secretions in the subglottic space.
- 2.No trauma to vocal cords when suctioning.
- 3.Convenient and safe.
- 4.Helping to reduce the rate of VAP (Ventilator-Associated Pneumonia) in the hospital.
- 5.Helping delay the early onset of VAP.
- 6.Size:5.5mm-8.5mm



Endobronchial Tube (Double Lumen)

Features

- 1.Pressure gauge valve series can be convenient for indicating the intracapsular.
- 2.pressure, airway injury to smallLow pressure high volume set of sac.
- 3.Imports of PVC material
- 4.Pure silica core one-way valve can be used in magnetic resonance imaging
- 5.Size:28FR/32FR/35FR/37FR/39FR/41FR.
- 6.Type:Left,Right.



Intubating Stylet (Bougie)

Features

- Made from medical grade PVC or Silicone.
- StyletSize:6FR/10Fr/14FR
- Bougie Size: 6FR/10Fr/14FR
- Type:Silod,Hollow



Endotracheal Tube Holder

Features

- 1.Fits various sizes of endotracheal tubes,from id 5 to id 10.
- 2.For single use only.



PVC Laryngeal Mask

Features

- 1.The tube is made of medical-grade PVC.
 - 2.The cuffs made of medical-grade silicone,and its specially designed shape coincides with the laryngopharynx well, reducing stimulation to patient body and imparoving the cuff seal.
 - 3.Suitable for adult, children and infant use.
- Size:1/1.5/2/2.5/3/4/5



Silicone Laryngeal Mask

Features

- 1.Made of medical-grade silicone.
 - 2.The cuff is made of medical-grade silicone,and its specially designed shape coincides with the laryngopharynx well, reducing stimulation to patient body and improving the cuff seal.
 - 3.Suitable for adult, children and infant use.
- 4.Size:1/1.5/2/2.5/3/4/5



Reusable Laryngeal Mask

Features

- 1.Made of medical-grade silicone.
- 2.Suitable for adult, children and infant use
- 3.Autoclave sterilization only;Can be reused up to 40 times.
- 4.Size: Size:1/1.5/2/2.5/3/4/5



Reinforced Laryngeal Mask

Features


- 1.Spiral reinforcement minimizes crushing or kinking.
- 2.Smooth, transparent and kink-resistant tube.
- 3.Suitable for adult, children and infant use.
- 4.Material: PVC /Silicone
- 5.Size: Size:1/1.5/2/2.5/3/4/5



Suction Catheter

Features

- 1.Made of non-toxic medical-grade PVC, clear and soft.
- 2.Color-coded connector for identification of different sizes.
- 3.Size:6FR-20FR



Guedel Airway

Features


- 1.Made form semirigid non-toxic PE color bide back.
- 2.Guedel style enclosed center channel.
- 3.Szie:40mm-120mm



Nasopharyngeal Airway

Features

- 1. PVC production, soft and comfortable.
- 2. For nasal ventilation only.
- 3. The tail end has a flared mouth type, flat mouth type, flat mouth thickened type and sliding ring type for selection.
- 4.Size:3.0mm-9.5mm



CHAPTER 3

URINARY

EXCELLENCE QUALITY INTEGRITY PROFESSIONAL

⌘ PRODUCTS



Urine Bag

Features

- 1.Made of medical grade PVC, PE, PP and ABS;
- 2.Consist of bag, connecting tube, taper connector, bottom outlet valva
- 3.Easy to read scale for quick determination of urine volume;
- 4.Four type available: without valve; screw valve; pushing valve; cross valve.



Urine Meter

Features

- 1.Made of medical grade PVC, PE, PP and ABS;
- 2.Consist of bag, connecting tube, taper connector, bottom outlet valva
- 3.Easy to read scale for quick determination of urine volume;
- 4.Four type available: without valve; screw valve; pushing valve; cross valve.



PVC Nelaton Catheter

Features

- 1.Medical grade PVC, transparent and soft
- 2.For urine drainage
- 3.Packed in sterilization pouch with ETO sterilization
- 4.Sizes available from Fr 8 - Fr 22



Latex Foley Catheter


Features

1. Medical grade natural rubber, soft and smooth
- 2.With silicone coating
- 3.2 way / 3 way available for options
- 4.Packed in sterilization pouch with ETO sterilization
5. Sizes available from Fr 8 - Fr 28

Silicone Foley Catheter

Features

- 1.Medical grade silicone,good bio-compatibility
- 2.Radiopaque line for x-ray visualization
- 3.2 way / 3 way available for options
- 4.Packed in sterilization pouch with ETO sterilization
5. Sizes available from Fr 8 - Fr 28



Urethral Catheter (Latex)

Features

- Material:Natural latex with silicone coated
- Size availabe from Fr6 - Fr32
- Length:400mm
- Type: one way



Condom Catheter (Latex)

Features

- Material: Natural latex
- Size: 20mm, 25mm, 30mm, 35mm, 40mm



Condom Catheter (Silicone)

Features

- Material: 100% silicone medical grade
- Size: 20mm, 25mm, 30mm, 35mm, 40mm



CHAPTER 4

- DIGESTIVE
- SURGICAL
- ICU PRODUCTS
- DRAINAGE SYSTEM

EXCELLENCE QUALITY
INTEGRITY PROFESSIONAL

DIGESTIVE PRODUCTS

Enema Bag

Features

Properties: Medical Materials & Accessories
Size: 2 L
Material: Other
Instrument classification: Class I
Product Name: enema bag kit
Material of bag: Rubber
Value: 2000 ml



Stomach Tube

Features


1. Medical grade PVC, transparent and smooth, silicone is optional
2. Radiopaque line for x-ray visualization
3. Packed in sterilization pouch or PE bag with ETO sterilization
4. Size available from Fr 12 - Fr 20
5. Length: 120 cm



Feeding Tube

Features

1. Medical grade PVC, transparent and smooth
2. Radiopaque line for x-ray visualization
3. Packed in sterilization or PE bag with ETO sterilization
4. Size available from Fr 5 - Fr 10 (Infant), Fr 12- Fr 20 (Adult)



Yankauer Handle And Suction Connecting Tube

Features

- 1.Handle made of K- Resin, tube made of medical grade PVC
- 2.Packed in sterilization pouch with ETO sterilization
- 3.4 types of handle available
- 4 handle type:crown tip with vent, crown tip without vent standard tip with vent, standard tip without vent
- 4. Customization is acceptable



Disposable Electrosurgical Pencil

Features

- Material:plastic
- Comfortable design and light weight
- Colour marked buttons,
- Hexagonal electrode connection,
- 3 meter cable, diameter 2.36 mm electrode
- 3 pin US-conncetor
- Blade can be changed to different length



Silicone Resuscitator

Features

- 1. Silicone resuscitator (except mask,oxygen tubing and reservior bag) can be autoclaved repeatedly at 134°C.
 - 2. It is equipped with a pressure limitation valve for Patient safety.
 - 3. Textured surface ensure a firm grip and provide effective ventilation.
 - 4. 100% latex free.
- Components: Mask, balloon with connector, oxygen tubing, reservoir bag, instruction sheet. All components placed in a PP box Or paper box.
- Size: Adult/L, Pediatric/M, Infant/S



PVC Resuscitator

Features

- 1. Due to being disposable,it can eliminate the possibility of cross-contamination. Any cleaning, disinfecting or sterilizing is not necessary for it.
 - 2. It is equipped with a pressure limitation valve for patient's safety.
 - 3. 100% latex free.
- Components: Mask, balloon with connector, oxygen tubing, reservoir bag, instruction sheet. All components placed in a PP box Or paper box.



CPR Mask

Features

- Material: PVC
- Size: Suitable for adult, children, infant
- Components: Mask, instruction sheet, PVC gloves and alcohol pad can be optional equipped.



First Aid Kit

Features

- Material: 600D oxford
- Can be configured according to customer requirements

Item	Specification
ICE bag	160g
PVC Gloves	M, powdered, nature
Non-woven Surgical Tape	1.25cmx5y, white, simple packing
Non-woven Surgical Tape	5cmx5y, white, simple packing
Wound Adhesive Plaster	72x19mm, PVC, brown color
PBT Elastic Bandage	3" x4m, bleach, cellophane packing
PBT Elastic Bandage	5cmx4m, cellophane packing
Gauze Bandage	5cmx5m
Gauze Bandage	7.5cmx5m

Item	Specification
Triangular Bandage	Mesh 46x30, 36"x36"x51"
Breathing Mask With One-Way Valve	20cmx20cm
Sterile Gauze Swab	Mesh 24x20, 10cmx 10cm-12p, 5pcs/bag
Sterile Gauze Swab	Mesh 24x20, 7.5cmx7.5cm-12p,5pcs/bag
Applicator	3", 6pcs/bag
Cotton Wool	3pcs/roll
Alcohol Prep Pad	65x30mm-2ply
Scissor	/
Tweezer	4#



PVC Wound Drainage System

Features

Material: PVC rubber, PP, PS, SS
Application: Used for abdominal, chest,breast and other parts of the fluid, pus and blood drainage.
Component: Container, drainage tube,connecting tube, needle, non-return valve.
Capacity: 200ml, 400ml, 500ml, 800ml
Size of Stainless Steel Needle: Fr7, Fr10, Fr12, Fr14, Fr15, Fr16, Fr18, Fr19
Packing: 1pc/Blister bag.



Silicone Wound Drainage System

Features

Material: silicone rubber, PP, PS, SS
Application: Used for abdominal, chest,breast and other parts of the fluid, pus and blood drainage.
Component: Container, drainage tube,connecting tube, needle, non-return valve.
Capacity: 200ml, 400ml, 500ml, 800ml
Size of Stainless Steel Needle: Fr7, Fr10, Fr12, Fr14, Fr15, Fr16, Fr18, Fr19.



MEDICAL SUPPLIES

EXCELLENCE QUALITY INTEGRITY PROFESSIONAL

88 PRODUCTS