



HANGZHOU WEHERE MEDICAL TECHNOLOGY CO., LTD.

#701, Block C, Wanda Plaza, Gongshu, Hangzhou, 310015, China

Tel.: +86-571-86822701 Cell.: +8613905717534

E-mail: moneyxd@weheremed.com

www.weheremedical.com

WEHERE  
MEDICAL  
*we here, we care*

# Gauze Series

Absorbent Gauze Roll  
Gauze Swabs  
Laparotomy Sponge  
Y/I Cut Gauze Swabs  
Gauze Cutting  
Gauze Ball  
Gauze Sheet  
Post-op Sponge  
Vaseline Gauze

Gauze Zigzag  
Calcium Alginate Wound Dressing (+AG)  
Honey Gauze Dressing



## Gauze Series

### Absorbent Gauze Roll




USP VII,BP 13,17,20 threads , white/green or blue gauze  
With/without X-ray detectable thread  
1-ply,2-ply or 4-ply

Jumbo roll 36"/48" x 1000/2000 m  
36"/48" x 5.5/11/55/100/110 yards

\*Other dimensions,width,length and package can be produced as per requirements

### Gauze Swabs



100% pure cotton USP VII,BP 13,17,20 threads , white,green or blue gauze  
With/without folded edge  
With/without X-ray detectable thread  
Available in STEAM,EO OR GAMMA sterile

5 x 5cm-8/12/16 ply  
7.5cm x 7.5cm-8/12/16 ply  
10cm x 10cm-8/12/16 ply  
10cm x 20cm-8/12/16 ply  
100/200 pcs/paper sleeve or box  
1/2/3/4/5/10 pcs/pouch or blister

\*Other dimensions,plies, and package can be produced as per requirements

### Laparotomy Sponge




BP 13,17,20 threads , white/green or blue gauze or non-woven  
With/without X-ray detectable thread or chip  
With/without metal or plastic ring,blue or white cotton loop  
Pre-washed or unwashed  
100pcs/polybag  
5pcs/blister with/without white or blue medical wrapping paper

|                     |                     |
|---------------------|---------------------|
| 20cm x 30cm-4/6 ply | 30cm x 30cm-4/6 ply |
| 45cm x 45cm-4/6 ply | 50cm x 60cm-4/6 ply |
| 90cm x 8cm-4ply     | 40cm x 40cm-4 ply   |
| 20cm x 20cm-12 ply  | 30cm x 30cm-12 ply  |
| 45cm x 45cm-12 ply  |                     |
| 18" x 18"-4ply      | 12" x 12"-4ply      |
| 4" x 18"-4ply       | 8" x 36"-4ply       |

\*Other dimensions,plies and package can be produced as per requirements

# Gauze Series

## Y/I Cut Gauze Swabs




100% pure cotton USP VII,BP 13,17,20 threads  
Y-cut or I-cut  
Sterile or non sterile  
100/200 pcs/sleeve or in paper-paper/paper-film/blister pack

100% cotton swab:  
10cm x 10cm-12/16 ply  
Non-woven swab:  
5cm/10cm-6 ply

\*Other dimensions ,plies and package can be produced as per requirements

## Gauze Cutting



100% pure cotton USP VII,BP 13,17,20 threads , white,green or blue gauze  
With/without X-ray detectable thread  
Sterile/Non sterile  
Packed in polybag/paper-paper/paper-film or blister pack

|             |             |             |
|-------------|-------------|-------------|
| 10cm x 10cm | 10cm x 20cm | 20cm x 20cm |
| 25cm x 25cm | 30cm x 30cm | 30cm x 40cm |

\*Other dimensions can be produced as per requirements

## Gauze Ball



100% pure cotton BP 17,20 threads  
With/without X-ray detectable thread  
In round,flat or hat shape  
Sterile or non sterile

Unfolded size:

|             |             |             |
|-------------|-------------|-------------|
| 10cm x 10cm | 15cm x 16cm | 20cm x 20cm |
| 25cm x 26cm | 30cm x 30cm | 30cm x 40cm |
| 40cm x 50cm |             |             |

\*Other dimensions and package can be produced as per requirements

# Gauze Series

## Gauze Sheet



100% pure cotton USP VII,BP 13,17,20 threads , white,green or blue gauze  
Sterile/Non sterile  
Packed in polybag/paper/paper-film bag or blister pack

|                  |                  |            |
|------------------|------------------|------------|
| 90cm x 18m       | 90cm x 27m       | 90cm x 36m |
| 90cm x 1m        | 90cm x 3m        | 90cm x 5m  |
| 90cm x 10m       |                  |            |
| 1 square meter   | 1/2 square meter |            |
| 1/4 square meter |                  |            |

\*Other dimension or package can be produced as per requirements

## Post-op Sponge



Covered with gauze/non woven  
Very high absorbency while with no lint  
Tissue in the middle,cotton inside or tissue inside only  
Machine/Manual folding  
Sterile/Non sterile  
Paper/paper-film pouch /blister pack

3" x 3"    3" x 4"    4" x 4"

\*Other dimension and package can be produced as per requirements

## Vaseline Gauze



Consists of gauze, paraffin/vaseline, oil emulsion  
21S/24X24  
200±20g/m2  
Available in different size  
1pc/individual aluminum foil bag packing;5pcs to 100pcs/mid box

5x5cm,10x10cm,10x20cm,10x40cm,  
10x700cm, 15x100cm, 15x200cm

\*Other specification can be produced as per requirements



# Gauze Series

## Gauze Zigzag



100%cotton medical absorbent surgical zigzag gauze is made of 100%cotton fiber by weaving and absorbent. It is white , soft and breath freely, with high absorbent capability , widely application in medical treatment and sanitation career.

Threads Mesh: 12x8 ,14X10,19X9 ,18x10,19x15, 24x20 , 25x17,26x18, 30x20 ,28x24, 30X22, 32X28,30x30 ,35X3040X30,40X50  
Length: 5m,15m,25m,50yards,50m, 100yards,100m,200m,500m

\*Other dimensions and package can be produced as per requirements

## Calcium Alginate Wound Dressing (+AG)



Alginate dressing is a kind of advanced medical wound dressing made of mixture of fiber and calcium ion from the natural seaweed. It can absorb 15-20 times their own weight in exudate.The calcium ions in the dressing interact with the sodium ions in the fluid of the wound. This reaction makes the fiber in the dressing to swell and partially dissolve into a gel.

Size: 5\*5/10\*10/2\*30 cm, 1 pc/pouch, 10-20 pcs/box

\*Other dimension ,weight or package can be produced as per requirements

## Honey Gauze Dressing



### Product performance

Manuka honey dressing is impregnated with 100% Manuka Honey from New Zealand with no additives, it is disposable sterile products, and consists of manuka honey and gauze. The dressing can protect the wound, effectively form and maintain the wet environment of the wound, and provide a suitable environment for wound healing.

| Dressing Size(cm) | Dressing Size(cm) |
|-------------------|-------------------|
| 5x5               | 10                |
| 10x10             | 10                |

Specifications can be customized

KING PHAR



# Non-woven Series

Non-woven Swabs

Drain Sponge

Non-woven Ball

Disposable Non-Woven shoe cover

FFP2 Protective Face Mask

Non-woven fabric roll

Non-Adherrent Pad

Isolation Gown/Coverall

Disposable PP Non woven Bed Cover

Medical Non-woven Wipes

Island Dressing

Fixing Type

Non-adherent foam dressing

Adult Diapers

Disposable Bath Towel

