



Hangzhou Wehere Medical Technology Co., Ltd

MAIN PRODUCTS CATALOGUE

主要产品目录

Hangzhou Wehere Medical Technology is a company on the medical health products manufacturing and trading in China. As a professional medical supplies nearly twenty years, our products have exported globally to Europe, Asia, America, etc., and help customers to succeed in their markets. We have three plants to produce different kinds of medical products. And we have set up longterm and close business relations with hundreds of domestic qualified manufactures.

On the basis of production management MDR and ISO13485 in our owner factories, guided by brand strategy, under the direct participation of American companies, we are committed to becoming an important member of the global medical and healthcare industry.

Add: #701, Block C, Wanda Plaza, Gongshu District, Hangzhou, China

Tel: 17815942351

Web: <https://www.weheremedical.com/>

E-mail: chen_chen@weheremed.com



1, oxygen mask with oxygen tube (氧气面罩)

Material: PVC

Mask Size: S/M/L/XL

Tube diameter: 5mm/6mm

Oxygen tube standard length: 2.1m





2, Tracheostomy mask (气切面罩)

Material: PVC

Mask Size: S/L/XL



3, Nebulizer mask （雾化面罩）

Material: PVC

Mask Size: S/M/L/XL

Jar size: 6cc, 8cc, 10cc, 20cc

Tube diameter: 5mm/6mm

Oxygen tube standard length: 2.1m

Color: green /transparent

Can be customized or OEM



4, Venturi mask with oxygen tube (文丘里面罩/氧浓度 6 色可调面罩)

Material: PVC

Mask Size: S/M/L/XL

Diluter: 6 oxygen concentration diluters

Tube diameter: 5mm/6mm

Oxygen tube standard length: 2.1m

Color: green /transparent

Can be customized or OEM



5, Venturi mask with oxygen tube (氧袋面罩)

Material: PVC

Mask Size: S/M/L/XL

Reservoir bag: 750ml / 1000ml

Tube diameter: 5mm/6mm

Oxygen tube standard length: 2.1m

Color: green /transparent

Can be customized or OEM



6, Mouthpiece nebulizer (T 型咬嘴雾化)

Material: PVC

Tube diameter: 5mm/6mm

Oxygen tube standard length: 2.1m

Color: green /transparent

Can be customized or OEM



7, Nasal oxygen cannula (鼻氧管)

Material: PVC

Tube diameter: 5mm/6mm

Oxygen tube standard length: 2.1m

Tube soft grade: common soft /super soft

Prong type: adult / pediatric / neonate

Color: green /transparent

Can be customized or OEM



8, single prong Nasal oxygen cannula (单腔吸氧管)

Material: PVC

Size; 8FR/ 10FR/ 12FR/14FR

Tube soft grade: common soft /super soft

Color: green /transparent

Can be customized or OEM



9, oxygen connecting tube (氧气连接管)

Material: PVC

Tube diameter: 5mm/6mm

Tube connector: standard / universal

Oxygen tube standard length: 2.1m

Color: green /transparent

Can be customized or OEM



10, Urinary catheter (导尿管)

Material: PVC

Size: Fr6-Fr18

Type: male / female

Tube Color: transparent

Can be customized or OEM



11, Urinary catheter bent at end (for special patient) (弯头导尿管)

Material: PVC

Size: Fr6-Fr18

Type: special patient

Tube Color: transparent

Can be customized or OEM



12, Suction catheter (吸痰管)

Material: PVC

Size: Fr5-Fr18

Tube Color: transparent

Can be customized or OEM



13, Feeding tube (喂食管)

Material: PVC

Size: Fr5-Fr18

Tube Color: transparent

Can be customized or OEM



14, Gastric tube (胃管)

Material: PVC

Size: Fr5-Fr18

Tube Color: transparent

Can be customized or OEM



15, Rectal tube (肛管)

Material: PVC

Size: Fr20-Fr32

Tube Color: transparent

Can be customized or OEM



16, Suction catheter kit (吸痰包)

Material: PVC

Size: Fr5-Fr18

Tube Color: transparent

Can be customized or OEM



17, Pediatric burette infusion set (小儿吊瓶输液器)

Material: PVC

Size: 150ml/100ml

Component is optional

Can be customized or OEM



18, Infusion set (输液器)

Material: PVC

Soft length: 1.2m-1.5m

Component is optional

Can be customized or OEM



19, Blood transfusion set (输血器)

Material: PVC

Soft length: 1.2m-1.5m

Component is optional

Can be customized or OEM



20, Light proof infusion set (避光输液器)

Material: PVC

Soft length: 1.2m-1.5m

Component is optional

Can be customized or OEM



21, Syringe (注射器)

Material: PVC

Size: 1ml/2ml/5ml/10ml/20ml/50ml/60ml

Can be customized or OEM

