

手术衣


SURGICAL GOWNS



Knitted cuff



Ultrasonic



Double sewing



Over-locking

"C&P" gowns has designed with safety, comfort and flexibility. We offer different fabrics which can meet EN13795 and AAMI Level3& Level4 performance requirements.


The variety of different material and models give you an option to choose the right gown for comfortable and protected procedure. Wide range of Surgical Gowns provides the advanced solution in different levels of protection to help protect patients and staff from the spread of infection during Surgery

Performance

- ⊙ Fluid Repellent
- ⊙ High Barrier performance
- ⊙ Cool and Breathable
- ⊙ Low Lint
- ⊙ Latex Free
- ⊙ Adjustable Velcro Neck Fasteners
- ⊙ Raglan Sleeves
- ⊙ Knitted Cuffs
- ⊙ Sterile gown is packed with 2 towels in a paper wrap and peelable bag


Standard Surgical Gown 标准手术衣

Size	Ref	L(cm)	W(cm)
M	GS010	115	150
L	GS011	125	157
XL	GS012	140	165
XXL	GS013	155	170




Reinforced Surgical Gown 加强型手术衣

Size	Ref	L(cm)	W(cm)
M	GR010	115	150
L	GR011	125	157
XL	GR012	140	165
XXL	GR013	155	170



T.U.R Surgical Gown 泌尿科手术衣

Size	Ref	L(cm)	W(cm)
XL	GT011	150	180



通用类手术


GENERAL SURGERY



General Pack



Universal Pack




Adhesive Drape



Fenestrated Drape



Laparotomy Drape




Fenestrated Drape

Universal Pack(REF:10101)

通用手术组合包

No.	REF	Description	QTY
1	8201	Back Table cover 140x190cm	1
2	8102	Mayo Stand cover 76x145cm	1
3	11402	Adhesive Side Drape 75x90cm	2
4		Adhesive Bottom Drape 150x180cm,	1
5		Adhesive Top Drape 150x240cm	1
6	90201	Hand Towel 30x40cm	4
7	90301	OP Tape Strip 10x50cm	1
8	GR011	Reinforced Surgical Gown, L size	2



Single Plain Drape

简易铺单


Size	REF
45x75cm	11001
75x90cm	11002
120x150cm	11003
60x60cm	11004



Adhesive Side Drape

贴胶边单


Size	Adhesive	Code
50x50cm	2.5X50cm	11401
75x90cm	4x90cm	11402
90x90cm	4x90cm	11403
100x150cm	4x100cm	11404
120x150cm	4x90cm	11405



Fenestrated Drape without Tape

洞巾不贴胶


Size	Fenestration	REF
50x50cm	Central Dia.7cm	11501
75x90cm	Central Oval 6x9cm	11502
90x90cm	Central Rectangular 10x12cm	11503
100x150cm	10x15cm Hole	11504
120x150cm	Central Square 10x10cm	11505



Fenestrated Drape with tape

洞巾贴膜


REF	Size(cm)	Description
11601	50x50	Central Dia.7cm,with tape
11602	75x90	Central Oval 6x9cm,with tape
11603	90x90	Central Rectangular 10x12cm,with tape
11604	100x150	Central Square 10x10cm,with tape
11605	120x150	Central Rectangular 12x15cm,with tape



Laparoscopy Pelviscopy Pack


(REF: 10201)腹腔镜手术包

No.	REF	Description	QTY
1	8201	Reinforced table cover140x190cm	1
2	8102	Mayo stand cover 76x145cm	1
3	11201	LaparascopyDrape,195x250x300cm,reinforcement 100x75cm,with 25x30cm fenestration,5cm adhesive around,Perineal fenestration20x10cm,with2.5cm tape,withflap over both fenestration	1
4	90301	Suture Bag	1
5	90201	Hand towel 30x40cm	4
6	GR011	Reinforced Surgical Gowns, L size	2



Laparoscopy Drape腹腔镜手术单


REF	Size(cm)	Description
11201	195x250x300	T shape,Absorbent area100x75cm,with 25x30cm fenestration,with 5cm adhesive ,Perineal fenestration20x10cm,with2.5cm tape,with flap over both fenestration
11202	200x260x310	T shape,25x30cm adhesive fenestration,2 instrument pouch,absorbent area100x75cm,4velcro tube holders



Laparotomy Pack (REF: 10301)

腹部手术包


No.	REF	Description	QTY
1	8202	Reinforced table cover150x220cm	1
2	8101	Mayo stand cover 58x137cm	1
3	11301	Laparotomy Drape,307x260x200cm,	1
4		Suture Bag	1
5		Utility Drapes with Adhesive 38 x 66 cm	4
6		1/2 Drapes 102 cm x 147 cm	2
7	90201	Hand towel 30x40cm	4
8	GR011	Reinforced Surgical Gowns, L size	2



Laparotomy Drape

腹部手术单


REF	Size(cm)	Description
11301	307x260x200	With armboard,Fabric reinforcement 75x75cm,with 10x30cm fenestration,2.5cm adhesive around,tube holders
11302	300x230	Fabric reinforcement 75x80cm with 12x25cm, fenestration, 2.5cm adhesive around



Thyroid (Struma) Drape

甲状腺手术单


REF	Size(cm)	Description
11701	190x300	Absorbent area 50x40cm, with diamond fenestration 12x12cm, 2.5cm adhesive around, 3 tube holder



Split Drape Pack(REF: 10401)


开叉手术包

No.	REF	Description	QTY
1	8201	Reinforced table cover140x190cm	1
2	8102	Mayo stand cover 76x145cm	1
3		Adhesive top drape,150x240cm,reinforcement 40x60cm	1
4		U Split drape,190x300cm,split 20x80cm,with 5cm adhesive and absorbent 60x120cm ,tube holder	1
5		Suture Bag	1
6	90201	Hand towel 30x40cm	4
7	GR011	Reinforced Surgical Gowns, L size	2



整形类手术

ORTHOPAEDICS



Hip Pack



Extremity Drape



Orthopaedic Pack



Hip Drape



Shoulder Arthroscopy Drape




Knee Arthroscopy Pack

Extremity Pack with gown(REF:20101)

四肢手术包


No.	Ref	Description	Qty
1	8201	Back Table Cover 140x190cm	1
2	8102	Mayo Stand Cover 76x145cm	1
3	21101	Extremity Drape,220x320cm	1
4		Utility Drape 75x75cm	1
5		OP Tape Strip 10x50cm	2
6	90201	Hand Towel 30x40cm	4
7		Impervious Split Drape , 140x175cm,10x50cm,U cut Adhesive Split	1
8	GR011	Reinforced Surgical Gown with 2pcs hand towel,XL	2



Shoulder Beach Chair Drape

肩关节镜手术单


REF	Size(cm)	Description
21201	260x400	reinforcement 75x50cm, Elastic fenestration 6cm Fluid collection pouch 75x90cm with 1 valve, 2 Velcro Tube holders



Shoulder Arthroscopy Drape

肩关节镜手术单


REF	Size(cm)	Description
21202	250x160	Split 15x50cm, Fluid collection pouch 27x115cm with 2 valve, Tape 5cm



Hand Pack with gown(REF:20301)

手部手术包


No.	Ref	Description	Qty
1	8201	Back Table Cover 140x190cm	1
2	8102	Mayo Stand Cover 76x145cm	1
3	21301	Hand Drape 270x300x160cm, Elastic fenestration 18x18cm with opening Dia.3cm, Absorbent reinforcement 50x80cm	1
4	11404	Medium Drape with Adhesive, 100x150cm	1
5	90301	OP Tape Strip 10x50cm	2
6	90201	Hand Towel 30x40cm	4
7	GR012	Reinforced Surgical Gown with 2pcs hand towel,XL	2



Hand Drapes

手部手术单


REF	Size(cm)	Description
21301	270x300x160	Elastic fenestration 18x18cm with opening Dia.3cm, Absorbent reinforcement 50x80cm, Gusseted hand table extension 160x90cm, Tube Holder



Hand &Foot Pack with

gown(REF:20401)手部足部手术包

No.	Ref	Description	Qty
1	8201	Back Table Cover 140x190cm	1
2	8102	Mayo Stand Cover 76x145cm	1
3	21401	Hand & Foot Drape 180x300cm, Elastic fenestration 18x18cm with opening Dia.3cm, Absorbent reinforcement 50x140cm	1
4	11404	Medium Drape with Adhesive, 100x150cm	1
5	90301	OP Tape Strip 10x50cm	2
6	90201	Hand Towel 30x40cm	4
7	GR012	Reinforced Surgical Gown with 2pcs hand towel,XL	2



Hand &Foot Drape

手部足部手术单


REF	Size(cm)	Description
21401	180x300	Elastic fenestration 18x18cm with opening Dia.3cm, Absorbent reinforcement 50x140cm



Knee Arthroscopy Pack with gown(REF:20501)

膝关节镜手术包

No.	Ref	Description	Qty
1	8201	Back Table Cover 140x190cm	1
2	8102	Mayo Stand Cover 76x145cm	1
3	21501	Arthroscopy T Drape 220x280x290cm, Fluid collection pouch 83x78cm, Circular Elastic fenestration 6cm,absorbent 75x120cm	1
4	90301	OP Tape Strip 10x50cm	2
5	90201	Hand Towel 30x40cm	4
6	90103	Impervious Stockinette,Medium 23x100cm	1
7		Suture Bag	1
8	GR012	Reinforced Surgical Gown with 2pcs hand towel,XL	2



Knee Arthroscopy Drapes with armboard cover膀胱镜手术单

REF	Size(cm)	Description
21501	220x280x290	Absorbent Reinforcement 75x120cm Circular Elastic fenestration Dia.6cm(Superior) Circular Elastic fenestration Dia.4cm(Inferior) Fluid collection pouch 83x78cm with 2 valve and flexwires Tube holders



Knee Arthroscopy Drape 膝关节镜手术单


REF	Size(cm)	Description
21502	220x320	Absorbent Reinforcement 75x120cm Circular Elastic fenestration Dia.6cm(Superior) Circular Elastic fenestration Dia.4cm(Inferior) Fluid collection pouch 60x115cm with 2 valve and flexwires 3 Tube holders



Hip Pack with gown(REF:20601)

臀部手术包


No.	Ref	Description	Qty
1	8201	Back Table Cover 140x190cm	1
2	8102	Mayo Stand Cover 76x145cm	1
3	21601	Hip Drape 250x290X320cm	1
4	90301	OP Tape Strip 10x50cm	2
5	90201	Hand Towel 30x40cm	4
6	90103	Impervious Stockinette,Medium 23x100cm	1
7		Suture Bag	1
8	GR012	Reinforced Surgical Gown with 2pcs hand towel,XL	2



Hip Drape With Pouch

臀部手术单


REF	Size(cm)	Description
21601	250x290x320	Armboard Cover,Absorbent Reinforcement 75x120cm Teardrop Elastic fenestration 20x17cm,Fluid collection pouch 60x120cm with 2 valve and flexwires Tube holders



Orthopaedic Pack with gown(REF:20701)

整形手术包


No.	Ref	Description	Qty
1	8201	Back Table Cover 140x190cm	1
2	8102	Mayo Stand Cover 76x145cm	1
3	21702	Ortho U Split drape 200x300cm	1
4	21701	Impervious Split Drape, 150x180cm,	1
5	90301	OP Tape Strip 10x50cm	1
6	90201	Hand Towel 30x40cm	4
7	90103	Impervious Stockinette,Medium 23x100cm	1
8		Suture Bag	1
9	GR012	Reinforced Surgical Gown with 2pcs hand towel,XL	2



Impervious U Drape

U形手术单


REF	Size(cm)	Description
21701	150x180	Split 10x53cm,with 5cm Adhesive



Ortho-U Drape

U形手术单


REF	Size(cm)	Description
21702	200x300	Split 10x100cm,with 5cm Adhesive,reinforcement 75x120cm,2tube holders



Vertical Isolation Drape

垂直隔离单

REF	Size(cm)	Description
21801	230x320	Incise flm 30x80cm, Fluid collection pouch 90x60cm with filter and valve, 2 Instrument pouch 25x28cm, Tape on top 3cm width, 2 Velcro Tube holder



妇产科、泌尿科手术


GYNAECOLOGY, OBSTETRICS & UROLOGY



T.U.R Drape



Under buttock Drape



C-Section Drape



Delivery Pack




Under buttock Drape

C-Section Drapes

剖腹产手术单


REF	Size(cm)	Description
31101	250x310x190	Full Incision Film 34x34cm, Fluid collection pouch with 2 suction port 75x90cm and forming wire, Absorbent zone 30x50cm



C-Section Drape

剖腹产手术单


REF	Size(cm)	Description
31102	250x310x190	Full Incision Film 34x34cm with Teardrop fenestration 25x20cm, Fluid collection pouch with 2 suction port 75x90cm and forming wire, Tube holders, Armboard Cover



Delivery Surgical Pack(REF:30201)

顺产手术包


No.	Ref	Description	Qty
1	8201	Back Table Cover 140x190cm	1
2	8102	Mayo Stand Cover 76x145cm	1
3	31203	UnderButtocks Drape, 102x112cm,with fluid collection pouch,absorbent zone	1
4	91101	Baby Blanket 60x90cm	2
5	90201	Hand Towel 30x40cm	4
6	11404	Utility Drape 100x150cm	1
7	90101	Leggings,75x120cm	2
8		Suture bag	1
9	GS012	Reinforced Surgical Gown with 2pcs hand towel,XL	2



Delivery Surgical Pack(REF:30202)

顺产手术包


No.	Ref	Name	Size	Qty
1		umbilical cord clamp	5cm L	2
2		shirt for woman in birth		1
3		stabling tracks	60x90cm	1
4		underwear		1
5		Sanitary napkin		1
6	GS011	Standard surgical gown,	L	2
7		Disposable bag	30x45cm	1
8		Facemask with earloop		1
9		bulb syringe		1
10		baby blanket	60x90cm	1
11		cotton swabs	10x10cm,8ply	5
12		Back table cover	150x100	1
13		Exam glove		2
14		ID bracelets		2



Under Buttock Drape

臀底单


REF	Size(cm)	Description
31201	90X120	Flim fabric reinforcement



Under Buttock Drape

臀底单


REF	Size(cm)	Description
31202	86X110	With Fluid Collection Pouch ,screen and suction port,forming wire, Absorbent zone 30x45cm



Under Buttock Drape

臀底单


REF	Size(cm)	Description
31203	102X112	Graduated Fluid collection pouch, screen and suction port,forming wire, Absorbent zone 30x45cm



Lithotomy Pack(REF:30301)


切石术手术包

No.	Ref	Description	Qty
1	8201	Back Table Cover 140x190cm	1
2	31301	Lithotomy T Drape, 150x220x75cm	1
3		Utility Drape with Adhesive 38 x 66cm	1
4	90101	Leggings75 cm x 120 cm	2
5	90201	Hand Towel 30x40cm	2



Lithotomy T Drape切石术手术单


REF	Size(cm)	Description
31301	150x200x75	Rectangular Adhesive Fenestration 10x15cm Absorbent Reinforcement 90x40cm,



T.U.R Pack(REF:30401)

泌尿科手术包


No.	Ref	Description	Qty
1	8201	Back Table Cover 140x190cm	1
2	31401	T.U.R drape 175x238x170cm	1
3	90505	Sterilization Wrap 100x100cm	1
4	90201	Hand Towel 30x40cm	1



T.U.R Drape

泌尿科手术单


REF	Size(cm)	Description
31401	175x238x170	Triangular fluid collection pouch,with screen, valve,with finger cot,fenestration Dia.7cm with tape zone of 12x18cm Tube holders



Cystoscopy Pack(REF:30501)

膀胱镜手术包

No.	Ref	Description	Qty
1	8201	Back Table Cover 140x190cm	1
2	31501	Cystoscopy T Drape 80x150x210	1
3	90505	Leggings 75x120cm	2
4	90201	Hand Towel 30x40cm	1



Cystoscopy T Drape

膀胱镜手术单

REF	Size(cm)	Description
31501	80x150x210	Perineal Fenestration 8x15cm, 22x22cm mesh screen, with tape, fabric reinforcement 75x120cm, 2 tube holders



眼科、耳鼻喉类手术

OPHTHALMOLOGY&ENT



Ophthalmic Pack



Eye Drape



ENT Split Drape



Eye Drape



Eye Drape



Eye Drape




Eye Drape

Cataract Pack (REF: 40101)

眼科手术包


No.	Ref	Name	Size	Qty
1	8205	Table Cover	150x140cm	1
2	41102	Ophthalmic Drape	100x120cm	1
3	GS011	Surgical Gown with towel	L	2
4		Syringe Luer slip	1ml, 27g	1
5		Syringe without needle Luer slip	3ml	2
6		Syringe without needle Luer slip	5ml	2
7		Syringe without needle Luer slip	10ml	1
8		Plastic bowl	60ml	1
9		Eye pad	5x7cm	1
10		Gauze Swab	7.5x7.5cm	10
11		Cellulose spears		5



Ophthalmic Pack (REF: 40102)


眼科手术包

No.	Ref	Name	Size	Qty
1	8205	Table Cover	150x140cm	1
2		Ophthalmic Drape	140x120cm	1
3	GS011	Surgical Gown with towel	L	2
4	90201	Hand towel	30x40cm	2
5		Syringe without needle Luer slip	2cc	1
6		Syringe without needle Luer slip	5cc	1
7		Syringe without needle Luer slip	10cc	1
8		Plastic container	150ml	1
10		Gauze Swab	10x10cm	5
11		Cellulose spears		5
12		Insulin syringe		1
13		Hydrodissection cannulas		1
14		Angled Reverse		1




Eye Drape眼科手术单

REF	Size(cm)	Description
41101	50x50	Round fenestration,with collection pouch
41102	100x120	Incise film 11x12cm,collection pouch25x30cm,fenestration size:9x9cm
41103	120x140	with one single fluid collection pouch,9x9cm fenestration,with incise,film
41104	150x200	with double fluid collection pouch,5x8cm oval fenestration,with incise, flexwire
41105	140x140	without fluid collection pouch,with oval fenestration6x9cm,incision film cut 7cm
41106	25x45cm	PE film incise, film,pouch



ENT Pack (REF:40201)耳鼻喉手术包

No	Ref	Name	Size	Qty
1	8205	Back Table Cover	150x140cm	1
2	8101	Mayo Stand Cover	58x137cm	1
3		Head Turban Drape	75x120cm	2
4	41201	ENT split drape.	190x300cm	2



ENT Split Drapes耳鼻喉手术单

REF	Size(cm)	Description
41201	190x300	with 8x70cm split, with adhesive, fabric reinforcement 75x120cm, 2 tube holder



神经类手术

NEUROLOGY








Cranitomy Drape



Cranitomy Drape


Neuro Pack(REF:50101)开颅手术包

No.	Ref	Name	Size	Qty
1	8201	Back Table Cover	190x140cm	1
2	8102	Mayo Stand Cover	76x145cm	1
3	51101	Craniotomy drape ,with incise film and pouch,oval fenestration	250x260 cm	1
4	90201	Absorbent Towel	30x40cm	2
5		Utility Drapes with adhesive tape	75X50 cm	2




Craniotomy Drape开颅手术单

REF	Size(cm)	Description
51101	250x260	50cm PE panel on both sides, Absorbent reinforced 50x160cm, Fluid Collection Pouch 50x76cm , Forming Wire, Mesh Strainer, Drainage Port, Oval Incise Film 26x20cm, Velcro Holding Tabs



Craniotomy T Drape开颅手术单

REF	Size(cm)	Description
51102	310x200x340	Absorbent reinforced 90x60cm, Fluid Collection Pouch 60x90cm , Forming Wire, Mesh Strainer, Drainage Port Oval Incise Film 26x20cm, Velcro Holding Tabs,Gusseted Anesthesia Panels



心脏心血管类手术

CARDIOVASCULAR & AUGIOGRAPHY




Augiography Drape

Angiography Drape Pack (REF60101)

心血管手术包

No.	Ref	Name	Size	Qty
1		Back Table Cover	190x140cm	1
2		Sterilization Wrap	90x90cm	1
3	61102	Femoral Angiography drape	220x320cm	1
4	90201	Hand Towel	30x40cm	4
5		Medium drape	75x90cm	1
6		Light cover round	Ø80cm	1
7		Light cover square,	100x150cm	1
8	GS012	Standard surgical gown	XL	2



Angiography Drape

心血管手术单

REF	Size(cm)	Description
61101	220x320	Dual Clear Side Panels 50cm width, Absorbent Reinforcement 80x100cm, Two Fenestration Ø14cm, two incision film with fenestration Ø 7.5cm



Radial Femoral Angiography Drape


心血管手术单

REF	Size(cm)	Description
61102	220x320	Dual Clear Side Panels 50cm width, Absorbent Reinforcement 100x75cm, Femoral fenestration Ø12cm, incision film with fenestration Ø 8cm, Radial fenestration oval 12cmx8cm, with fenestration Ø 6x4cm




Femoral Angiography Drape with pouch 心血管手术单

REF	Size(cm)	Description
61103	240x370	Dual Clear Side Panels 50cm width, Absorbent Reinforcement 125x130cm, 2 Collection pouch 125x15cm, Femoral fenestration Ø12cm, incision film with fenestration Ø6cm, Radial fenestration 12cm, with fenestration Ø9cm



Brachial Angiography Drape 心血管手术单

REF	Size(cm)	Description
61104	80x110	Absorbent Reinforcement 30x60cm, Round central Fenestration of 7cm with incise film




Cardiovascular Pack 心脏手术包

No.	Ref	Name	Size	Qty
1	8102	Back Table Cover	190x140cm	1
2	90505	Sterilization Wrap	100x100cm	1
3	61201	Cardiovascular Drape	200x300x350cm	1
4	90201	Hand Towel	30x40cm	4
5		Perineal Drape	190x240 cm	1
6		Foot Cover		2
7		Utility Drape with 75x75cm tape		2
8	90301	OP tape	10x50cm	1
9		Suture Bag		1
10	GS012	Standard surgical gown	XL	2

Cardiovascular Drape 心脏手术单

REF	Size(cm)	Description
61201	200x300x350	Chest Incise Film 33x38cm, Leg Incise Film 55x96cm, 2 Fluid Collection Pouch, Absorbent Reinforcement 110x210cm, Velcro Holding Tabs, Armboard Cover



牙科手术包

DENTAL PACK

NO.	REF	Description	QTY
1	11402	Drape 90x90cm with lateral adhesive	1
2	11002	Drapes 75x90 cm	2
3	11504	Drape 100x150 cm, fenestration 10 cm	1
4		Protections for tubing 120x8 cm	2
5		Drape transparent 30x30 cm, with adhesive strip	1
6		Nonwoven swaps 5*5cm(Pack of 10 swap)	10
7		Transparent Plastic bag	1
8	90301	Adhesive strip 50x5 cm	2
9		PE Lamp cover	2
10		Plastic cover for radio tube (brilliance amplifactor)	1
11	90203	Hand towel	4
12	GS011	Surgical Gowns SMS size L	2
13	0601	Isolation gown with elastic cuff 23 g	1
14	0102	Surgeon caps - with tie	2
15	0101	Nurse caps with elastic	2
16	0202	Surgical face masks with tie	2
17	0301	Shoe covers (PE)	6



设备仪器罩


EQUIPMENT COVER



Mayo stand Cover



Mayo stand Cover



Back Table Cover



Banded Bag