



WEHERE MEDICAL

Wound Care Specialist

contents

Product Name	Page
Transparent Dressing	1
Incise Drape	3
Transparent Film Fixation Roll	4
Non Woven Fixation Roll	5
Transparent Wound Dressing	7
Adhesive Wound Dressing	9
Silver Antibacterial Dressing	11
Alginate Dressing	12
Hydrocolloid Dressing	15
Foam Dressing	17
Silicone Wound Contact Dressing	19
Silicone Tape	21
Cannular Fixation Device	23
NPWT Dressing Kit	25
Jumbo Rolls	26
Product Size	28

Transparen Dressing




Customized
Service
Available



Transparent IV dressing is designed to provide secure and reliable protection for IV sites, allowing for quick and easy access to the site. It is made from soft and flexible PU film that is transparent and breathable, providing a comfortable and secure fit. It is also waterproof and bacteria barrier, ensuring that the IV site stays dry and protected.



Incise Drape



The surgical incise drape provides optimal protection and reliable visibility during surgical procedures. It is permeable to moisture vapour, preventing fluids or debris from entering the incision site.


Contain Iodine and could keep the operation part sterile during operation.

Size(cm)	Packing
10x14	20pcs/box , 40boxes/ctn
12x25	20pcs/box , 20boxes/ctn
15x28	10pcs/box , 20boxes/ctn
20x30	20pcs/box , 20boxes/ctn
28x30	20pcs/box , 20boxes/ctn
30x45	20pcs/box , 10boxes/ctn
38x41	10pcs/box , 20boxes/ctn
45x45	20pcs/box , 10boxes/ctn
60x45	10pcs/box , 20box/ctn
90x60	10pcs/box , 20boxes/ctn
90x45	10pcs/box , 10boxes/ctn
90x45	10pcs/box , 10boxes/ctn
Shape, size and packaging can be customized.	



Transparent Film Fixation Roll

- Allow shower and swim
- Barrier to outside contamination
- Could be cut as desired
- Hypoallergenic
- Non-irritating
- Secure fixation
- Breathable
- High elasticity



Tattoo Aftercare

Non Woven Fixation Roll

- Hypoallergenic
- Non-irritating
- Secure fixation
- Could be cut as desired
- Comfortable to wear



Size(cm)	Packing
5cmx10m	1roll/box , 40boxes/ctn
7cmx10m	1roll/box , 30boxes/ctn
10cmx10m	1roll/box , 30boxes/ctn
15cmx10m	1roll/box , 20boxes/ctn
20cmx10m	1roll/box , 20boxes/ctn

Shape, size and packaging can be customized.





- Breathable
- Transparent
- Waterproof
 - Bacteria barrier
 - Stretchable
 - Hypoallergenic
 - Non irritating
 - Diversified sizes
 - High absorbent pad
 - Non-adherent to wound



Could combine with different core, such as high absorbent pad, alginate, hydrocolloid, foam and so on.



Could be made into various sizes and shapes.


Transparent Wound Dressing

Transparent Wound Dressing is made of a thin, transparent PU film plus a high absorbent pad that allows for easy monitoring of the wound site. It is highly breathable and conforms to the body for maximum comfort and protection. It is waterproof, bacteria-resistant, and provides a secure seal to keep the wound clean and protected. The absorbent pad could absorb wound exudate effectively and won't adhere to wound.





Adhesive Wound Dressing



Could be made into different sizes and shapes.


- Skin friendly
- High absorbent
- Non-adherent to wound
- Breathable
- Comfortable to wear




The adhesive wound dressing is made with a comfortable and breathable spunlace non woven plus a high absorbent pad that is gentle on the skin. It is designed to provide a secure, breathable barrier that helps to keep the wound clean and promote healing. The dressing is highly absorbent and adheres securely to the skin, making it ideal for a wide range of wound types. The absorbent pad could absorb wound exudate effectively and won't adhere to wound.



Silver Antibacterial Dressing is a revolutionary new product made with silver ion technology. It is designed to provide long-lasting protection against bacteria and fungi, while also promoting healing and preventing infection. It is high absorbent and non-adherent to wound site.



Silver Antibacterial Dressing



Alginate Dressing




Silver Alginate Dressing



Alginate dressing is an advanced wound care product designed to help promote wound healing. It is made from seaweed and is a natural, absorbent material that helps to draw out infection and reduce pain and discomfort. This dressing is easy to use and can be cut to fit for various wound sizes.



- High absorbent pad
- Non-irritating
- Good gel forming ability
- Non-adherent to wound
- Diversified sizes



- Maintain a moist wound environment
- Non-adherent to wound
- Diversified shape and size.

Hydrocolloid dressing is used to absorb wound exudate, promote moist wound healing, and protect the wound from further damage. The dressing is easy to apply and



remove, and can be used on a variety of wounds, including abrasions, lacerations, and ulcers. The hydrocolloid particles also help to reduce pain and discomfort, and promote faster healing.

Hydrocolloid Dressing





Foam Dressing



Ag



Foam Dressing is designed to provide superior protection and comfort for wound site. Its unique foam technology helps to create a breathable barrier that locks in wound exudation and keeps out dirt and bacteria. The foam dressing is highly

absorbent, allowing it to quickly and effectively draw away exudate and help to reduce the risk of infection. With silicone adhesive, it stays well on skin, while non irritating to skin and won't hurt skin when remove.

- High absorbent
- Breathable
- Hypoallergenic
- Non adherent to wound
- Non stick hair
- Waterproof
- Bacteria Barrier
- Diversified shape and size
- Comfortable

Silicone Wound Contact Dressing



Size(cm)	Packing	
5cmx7cm	1pc/pouch	10pcs/box
7.5cmx10cm	1pc/pouch	10pcs/box
10cmx10cm	1pc/pouch	10pcs/box
10cmx15cm	1pc/pouch	10pcs/box
20cmx30cm	1pc/pouch	10pcs/box

Size and packing could be customized.



- Painless Removal
- Gentle to Fragile Skin
- Non Toxic & Non Allergy
- Allow Exudate Pass Through Freely
- Reduce Damage To Wounds When Changing Dressings


Silicone Tape

Surgical silicone tape is designed for long-term wear and can be used to secure bandages, dressings, and other medical devices. The tape is hypoallergenic, non-toxic, and safe to use on sensitive and fragile skin. It can be used repeatedly without losing its stickiness.

- Hypoallergenic
- Comfortable
- Non adherent to wound
- Prevent scar formation
- Reusable
- Non stick hair
- Waterproof
- Could be cut as desired



Size(cm)	Packing	
2cmx150cm	1pc/box	100boxes/ctn
2cmx300cm	1pc/box	100boxes/ctn
4cmx150cm	1pc/box	100boxes/ctn
5cmx150cm	1pc/box	100boxes/ctn
Size and packaging could be customized.		



3.5cm x 9cm



For thin drainage tube, PICC,
CVC tube and indwelling needle fixation



Epidural anesthesia catheter fixation,
indwelling needle fixation



Central venous catheter fixation



Central venous catheter fixation

Heparin cap fixation of indwelling needle

2cm x 7.3cm

Fixation of gastric tube and nasogastric feeding tube



6.5cmx10.5cm

For fixing PICC and CVC catheter, central venous catheter fixation



6.5cmx10.5cm

For various drainage tubes and catheter fixation



3cmx9.2cm



8cm x 11cm



3.5cm x 9cm



Cannular Fixation Device



6cm x 8.5cm

nasogastric tubes,gastric tubes and oxygen tubes (both legs)



6.5cm x 9.5cm

nasogastric tubes,gastric tubes and oxygen tubes (single leg)

Product Features

Firm fixation, avoid pipeline slippage,
reduce patient pain, replace traditional
tape fixation and make clinical
operation more convenient.

6cm x 7cm

Fixation of indwelling needle




Negative Pressure Wound Therapy Dressing Kits.

Black Negative Foam	1PC	10x7.5cm, 10x15cm, 15x20cm, 15x25cm thickness : 2cm/2.5cm/3cm
Transparent Incise Film	1PC	15cmx20cm , 20cmx30cm , 25cmx30cm , 30cmx30cm , 30x45cm
Suction Interface with Tubing	1PC	0.5m , 0.8m , 1.5m

Shape size and packaging could be customized.



NPWT Dressing Kit



outer diameter : 10cm

Negative Pressure Wound Therapy (NPWT) is used for treating chronic and difficult-to-heal wounds. A key component of NPWT is the wound dressing material, which have to

maintain a seal with the skin and effectively transmit negative pressure to the wound bed. We could produce as your demand to suit different market needs.

Jumbo Rolls



Size(cm)	Packing
30cmx100m	2rolls/ctn
45cmx200m	1roll/ctn
100cmx300m	1roll/ctn
Shape, size and packaging can be customized.	

We could not only provide final products, semi-finished products, but also jumbo rolls according to your special requirements to meet different market demands.



Transparent Dressing

Size(cm)	Packing
4cmx5cm	1pc/pouch,50pcs/box
6cmx7cm	1pc/pouch,50pcs/box
10cmx12cm	1pc/pouch,50pcs/box
12cmx12cm	1pc/pouch,50pcs/box
10cmx15cm	1pc/pouch,20pcs/box
12cmx25cm	1pc/pouch,20pcs/box
15cmx20cm	1pc/pouch,20pcs/box
20cmx30cm	1pc/pouch,20pcs/box
10cmx20cm	1pc/pouch,20pcs/box
5cmx5.7cm	1pc/pouch,50pcs/box
7cmx9cm	1pc/pouch,50pcs/box
8.5cmx10.5cm	1pc/pouch,50pcs/box
7cmx10.5cm	1pc/pouch,50pcs/box
9cmx12cm	1pc/pouch,50pcs/box
10cmx11.5cm	1pc/pouch,50pcs/box

Hydrocolloid Dressing

Silicone Foam Dressing Border / Without Border

Silver Foam Wound Dressing

Foam Dressing / Sliver Foam Dressing

Size(cm)	Packing
5cmx5cm	1pc/pouch,20pcs/box
7.5cmx7.5cm	1pc/pouch,20pcs/box
10cmx10cm	1pc/pouch,20pcs/box
12.5cmx12.5cm	1pc/pouch,10pcs/box
15cmx15cm	1pc/pouch,10pcs/box
14cmx18cm (sacrum)	1pc/pouch,10pcs/box
15cmx18cm (sacrum)	1pc/pouch,10pcs/box
17cmx17.5cm (sacrum)	1pc/pouch,10pcs/box
23cmx23cm (sacrum)	1pc/pouch,10pcs/box

Transparent Wound Dressing

Adhesive Non-woven Wound Dressing

Alginate Dressing/Silver Alginate Dressing

Silver Antibacterial Dressing

Size(cm)	Packing
5cmx5cm	1pc/pouch,50pcs/box
6cmx7cm	1pc/pouch,50pcs/box
9cmx10cm	1pc/pouch,20pcs/box
9cmx15cm	1pc/pouch,20pcs/box
10cmx10cm	1pc/pouch,20pcs/box
10cmx12cm	1pc/pouch,20pcs/box
10cmx15cm	1pc/pouch,20pcs/box
10cmx20cm	1pc/pouch,20pcs/box
10cmx30cm	1pc/pouch,20pcs/box
20cmx35cm	1pc/pouch,20pcs/box
12cmx12cm(Drainage)	1pc/pouch,20pcs/box
12cmx14cm(Drainage)	1pc/pouch,20pcs/box
6.3cmx9.5cm(Eye)	1pc/pouch,50pcs/box
6cmx8cm(Eye)	1pc/pouch,50pcs/box
7cmx9cm(Eye)	1pc/pouch,50pcs/box

We could produce as your demand if
above size do not fulfill your requirements.